

Datasheet for 209-401-C24**IL-29 Antibody****Overview**

Description:	Anti-Human IL-29 (RABBIT) Antibody - 209-401-C24
Item No.:	209-401-C24
Size:	100 µg
Applications:	WB
Reactivity:	Human
Host Species:	Rabbit

Product Details

Background:	IL-29 (Interleukin-29) is a secreted cytokine with immunomodulatory activity. It may play a role in antiviral immunity. It up-regulates MHC class I antigen expression. IL-29 is a ligand for the heterodimeric class II cytokine receptor composed of IL10RB and IL28RA. The ligand/receptor complex seems to signal through the Jak-STAT pathway. Its induction is by viral infections or double stranded RNA. Anti-IL-29 antibody is ideal for investigators involved in Immunology research.
Synonyms:	rabbit anti-IL-29 antibody, rabbit anti-Interferon lambda-1 antibody, rabbit anti-IFN-lambda-1 antibody, Interferon lambda-1, Cytokine Zcyto21, Interleukin-29, IL-29, IFNL1, IL29, ZCYTO21
Host Species:	Rabbit
Clonality:	Polyclonal
Format:	IgG

Target Details

Gene Name:	IFNL1
Reactivity:	Human
Immunogen Type:	Recombinant Protein
Immunogen:	This IgG fraction antibody was prepared from rabbit antiserum after repeated immunizations with mature length recombinant human IL-29 protein produced in E.coli.

Purity/Specificity: Anti-IL-29 is an IgG fraction antibody purified from monospecific antiserum by a multi-step process which includes delipidation, salt fractionation and ion exchange chromatography followed by extensive dialysis against the buffer stated above. This antibody is specific for human IL-29 protein. A BLAST analysis was used to suggest cross-reactivity with IL-29 from human sources based on 100% homology with the immunizing sequence. Based on high to 100% homology, there is a chance of cross-reactivity to IL-29 from a wide variety of animals. Cross-reactivity with IL-29 from other sources has not been determined.

Relevant Links:

- [UniProtKB - Q8IU54](#)
- [NCBI - NP_742152.1](#)
- [GeneID - 282618](#)

Application Details

Tested Applications:	WB
Application Note:	This purified antibody has been tested in western blotting and suitable for ELISA. Specific conditions for reactivity should be optimized by the end user. Expect a band approximately 19.9 kDa in size corresponding to the mature human IL-29 protein by western blotting in appropriate cell lysate or extract.
Assay Dilutions:	All assays should be optimized by the user. Recommended dilutions (if any) may be listed below.
ELISA:	1:10,000
WB:	1:1000

Formulation


Physical State:	Lyophilized
Concentration:	1.0 mg/mL by UV absorbance at 280 nm
Buffer:	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2
Preservative:	None
Stabilizer:	None
Reconstitution Volume:	100 µL
Reconstitution Buffer:	Restore with deionized water (or equivalent)

Shipping & Handling

Shipping Condition: Ambient

Storage Condition:	Store anti-IL-29 antibody at 4° C prior to restoration. For extended storage aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use.
Expiration:	Expiration date is six (6) months from date of receipt.

Images



Western Blot

Anti-human IL-29 antibody in western blot shows detection of recombinant human IL-29 raised in E.coli. Recombinant protein (0.1 µg, 19.9 kDa) was loaded onto and resolved by SDS-PAGE, then transferred to nitrocellulose. The membrane was blocked with 1% BSA in TBST for 30 min at RT, followed by incubation with Rockland's, Inc. Anti-Human IL-29. After washing, membrane was probed with secondary antibody Dylight™ 649 Conjugated Anti-Rabbit IgG (H&L) (Goat) Antibody (611-143-122) diluted 1:20,000 in blocking buffer (p/n MB-070) for 30 min. at RT. Data was collected using Bio-Rad VersaDoc® 4000 MP imaging system.

Disclaimer

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 5199, Limerick, Pennsylvania, USA.