

Datasheet for 209-401-C27**BDNF Antibody****Overview**

Description:	Anti-BDNF (RABBIT) Antibody - 209-401-C27
Item No.:	209-401-C27
Size:	100 µg
Applications:	WB
Reactivity:	Human
Host Species:	Rabbit

Product Details

Background:	BDNF is a member of the nerve growth factor family of trophic factors. In the brain BDNF has a trophic action on retinal, cholinergic, and dopaminergic neurons, and in the peripheral nervous system it acts on both motor and sensory neurons. Some protein domains of BDNF are identical with those of NGF and another neurotrophic factor, designated NT-3 (neurotrophin-3). It exists as monomers and homodimers, and binds to NTRK2/TRKB. Polyclonal antibodies raised against murine NGF have been shown to cross-react with both NT-3 and BDNF. The propeptide is N-glycosylated and glycosulfated. BDNF is converted into mature BDNF by plasmin (PLG). It is expressed in brain, and highly expressed in hippocampus, amygdala, cerebral cortex and cerebellum; it is also expressed in heart, lung, skeletal muscle, testis, prostate and placenta.
Synonyms:	rabbit anti-BDNF, rabbit anti-brain derived neurotropic factor, Brain-derived neurotrophic factor, Abreineurin cytokine
Host Species:	Rabbit
Clonality:	Polyclonal
Format:	IgG

Target Details

Gene Name:	BDNF
Reactivity:	Human
Immunogen Type:	Recombinant Protein

Immunogen:	This IgG fraction antibody was prepared from rabbit antiserum after repeated immunizations with recombinant truncated human BDNF protein produced in E.coli.
Purity/Specificity:	This product is an IgG fraction antibody purified from monospecific antiserum by a multi-step process which includes delipidation, salt fractionation and ion exchange chromatography followed by extensive dialysis against the buffer stated above. This antibody is specific for human BDNF protein. A BLAST analysis was used to suggest cross-reactivity with BDNF from human sources based on 100% homology with the immunizing sequence. Based on high to 100% homology, there is a chance of cross-reactivity to BDNF from a wide variety of animals. Cross-reactivity with BDNF from other sources has not been determined.
Relevant Links:	• UniProtKB - P23560 • NCBI - CAA62632.1 • GenelD - 627

Application Details

Tested Applications:	WB
Application Note:	This purified antibody has been tested in western blotting and suitable for ELISA. Specific conditions for reactivity should be optimized by the end user. Expect a band approximately 27 kDa in size corresponding to the mature human BDNF protein by western blotting in appropriate cell lysate or extract.
Assay Dilutions:	All assays should be optimized by the user. Recommended dilutions (if any) may be listed below.
ELISA:	1:10,000
WB:	1:1000


Formulation

Physical State:	Lyophilized
Concentration:	1.0 mg/mL by UV absorbance at 280 nm
Buffer:	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2
Preservative:	None
Stabilizer:	None
Reconstitution Volume:	100 µL
Reconstitution Buffer:	Restore with deionized water (or equivalent)

Shipping & Handling

Shipping Condition:	Ambient
Storage Condition:	Store vial at 4° C prior to restoration. For extended storage aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use.
Expiration:	Expiration date is six (6) months from date of receipt.

Images



Western Blot

Anti-human BDNF antibody in western blot shows detection of recombinant human BDNF raised in E.coli. Recombinant truncated (0.1 µg, 27 kDa) protein was loaded onto and resolved by SDS-PAGE, then transferred to nitrocellulose. The membrane was blocked with 1% BSA in TBST for 30 min at RT, followed by incubation with Rockland's, Inc. Anti-Human BDNF. After washing, membrane was probed with secondary antibody Dylight™ 649 Conjugated Anti-Rabbit IgG (H&L) (Goat) Antibody (611-143-122) diluted 1:20,000 in blocking buffer (p/n MB-070) for 30 min. at RT. Data was collected using Bio-Rad VersaDoc® 4000 MP imaging system.

Disclaimer

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 5199, Limerick, Pennsylvania, USA.