

**Datasheet for 209-403-B96****IL-9 Antibody Peroxidase Conjugated****Overview**

|                      |                                                                       |
|----------------------|-----------------------------------------------------------------------|
| <b>Description:</b>  | Anti-Human IL-9 (RABBIT) Antibody Peroxidase Conjugated - 209-403-B96 |
| <b>Item No.:</b>     | 209-403-B96                                                           |
| <b>Size:</b>         | 100 µg                                                                |
| <b>Applications:</b> | WB                                                                    |
| <b>Reactivity:</b>   | Human                                                                 |
| <b>Host Species:</b> | Rabbit                                                                |

**Product Details**

|                      |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
|----------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Background:</b>   | IL-9 is a cytokine that acts as a regulator of a variety of hematopoietic cells. This cytokine stimulates cell proliferation and prevents apoptosis. It functions through the interleukin 9 receptor (IL9R), which activates different signal transducer and activator (STAT) proteins and thus connects this cytokine to various biological processes. The gene encoding this cytokine has been identified as a candidate gene for asthma. Genetic studies on a mouse model of asthma demonstrated that this cytokine is a determining factor in the pathogenesis of bronchial hyperresponsiveness. |
| <b>Synonyms:</b>     | rabbit anti-IL-9 peroxidase conjugated antibody, rabbit anti-interleukin-9 HRP conjugated antibody, Interleukin 9, IL 9, IL9, T cell growth factor 3, T cell growth factor p40, TCGF 3, Mast cell growth factor, MCGF, Megakaryoblast growth factor, p40 cytokine                                                                                                                                                                                                                                                                                                                                    |
| <b>Host Species:</b> | Rabbit                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| <b>Conjugate:</b>    | Peroxidase (HRP)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>Clonality:</b>    | Polyclonal                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Format:</b>       | IgG                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |

**Target Details**

|                        |                     |
|------------------------|---------------------|
| <b>Gene Name:</b>      | IL9                 |
| <b>Reactivity:</b>     | Human               |
| <b>Immunogen Type:</b> | Recombinant Protein |

|                            |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
|----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Immunogen:</b>          | This purified antibody was prepared from whole rabbit serum produced by repeated immunizations with full length recombinant human IL-9 protein.                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>Purity/Specificity:</b> | This product is an IgG fraction antibody purified from monospecific antiserum by a multi-step process which includes delipidation, salt fractionation and ion exchange chromatography followed by extensive dialysis against the buffer stated above. This purified antibody has been heated to 56°C for 30 minutes. In ELISA and other immunoreactive assays, this antibody will recognize both native and recombinant human IL-9 in cell supernatants and certain body fluids. A control of similarly diluted normal rabbit IgG is recommended. |
| <b>Relevant Links:</b>     | <ul style="list-style-type: none"><li>• <a href="#">UniProtKB - P15248</a></li><li>• <a href="#">NCBI - NP_000581.1</a></li><li>• <a href="#">GenelD - 3578</a></li></ul>                                                                                                                                                                                                                                                                                                                                                                         |

## Application Details

|                             |                                                                                                                                                                                                                                                                                                           |
|-----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Tested Applications:</b> | WB                                                                                                                                                                                                                                                                                                        |
| <b>Application Note:</b>    | This purified antibody has been tested in western blotting and suitable for ELISA. By western blot a band approximately 15 kDa in size corresponding to human IL-9 protein is expected in the appropriate cell lysate or extract. Specific conditions for reactivity should be optimized by the end user. |
| <b>Assay Dilutions:</b>     | All assays should be optimized by the user. Recommended dilutions (if any) may be listed below.                                                                                                                                                                                                           |
| <b>ELISA:</b>               | 1:10,000-1:50,000                                                                                                                                                                                                                                                                                         |
| <b>IHC:</b>                 | 1:500-1:2,500                                                                                                                                                                                                                                                                                             |
| <b>WB:</b>                  | 1:1000-1:5,000                                                                                                                                                                                                                                                                                            |


## Formulation

|                               |                                                                        |
|-------------------------------|------------------------------------------------------------------------|
| <b>Physical State:</b>        | Lyophilized                                                            |
| <b>Concentration:</b>         | 1.0 mg/mL                                                              |
| <b>Buffer:</b>                | 0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2             |
| <b>Preservative:</b>          | 0.01% (w/v) Gentamicin Sulfate. Do NOT add Sodium Azide!               |
| <b>Stabilizer:</b>            | 10 mg/mL Bovine Serum Albumin (BSA) - Immunoglobulin and Protease free |
| <b>Reconstitution Volume:</b> | 100 µL                                                                 |
| <b>Reconstitution Buffer:</b> | Restore with deionized water (or equivalent)                           |

## Shipping & Handling

|                            |                                                                                                                                                                                                                                                                                                                                                 |
|----------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Shipping Condition:</b> | Ambient                                                                                                                                                                                                                                                                                                                                         |
| <b>Storage Condition:</b>  | Store vial at 4° C prior to restoration. For extended storage aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. |
| <b>Expiration:</b>         | Expiration date is one (1) year from date of receipt.                                                                                                                                                                                                                                                                                           |

## Images



### Western Blot

Western Blot of Rabbit anti-IL-9 Antibody Peroxidase Conjugated. Lane 1: Human IL-9. Load: 50 ng per lane. Secondary antibody: IL-9 Antibody Peroxidase Conjugated at 1:1,000 for 30 min at RT. Block: MB-070 for 30 min RT. Predicted/Observed size: 14 kDa, 14 kDa.

## Disclaimer

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 5199, Limerick, Pennsylvania, USA.