

Datasheet for 210-301-B89**Esrp-1 Antibody****Overview**

Description:	Anti-Esrp-1 (MOUSE) Monoclonal Antibody - 210-301-B89
Item No.:	210-301-B89
Size:	100 µg
Applications:	WB
Reactivity:	Mouse
Host Species:	Mouse

Product Details

Background:	Epithelial splicing regulatory protein 1 (Esrp-1) is an mRNA splicing factor that regulates the formation of epithelial cell-specific isoforms. It specifically regulates the expression of FGFR2-IIIb, an epithelial cell-specific isoform of FGFR2, and also regulates the splicing of CD44, CTNND1, ENAH, 3 transcripts that undergo changes in splicing during the epithelial-to-mesenchymal transition (EMT). EsrpP-1 acts by directly binding specific sequences in mRNAs. It binds the GU-rich sequence motifs in the ISE/ISS-3, a cis-element regulatory region present in the mRNA of FGFR2.
Synonyms:	mouse anti-Esrp-1 antibody, mouse anti-Esrp1 antibody, Rbm35a, Epithelial splicing regulatory protein 1, RNA-binding protein 35A, RNA-binding motif protein 35A
Host Species:	Mouse
Clonality:	Monoclonal
Clone ID:	27H12.F3.H5
Format:	IgG2b

Target Details

Gene Name:	Esrp1
Reactivity:	Mouse
Immunogen Type:	Recombinant Protein

Immunogen:	Anti-Esrp-1 was produced by repeated immunizations of full length recombinant mouse Esrp-1 fusion protein.
Purity/Specificity:	This product is an IgG fraction antibody purified by Protein-A chromatography followed by extensive dialysis against the buffer stated above. This antibody reacts with mouse Esrp-1 protein. A BLAST analysis of the immunizing protein sequence shows 100% homology with Esrp-1 from mouse and a 91% sequence homology with Esrp-1 from human, pig, rat, opossum, horse, cattle, panda, dog, and chimpanzee. The binding epitope of this monoclonal antibody has not been mapped.
Relevant Links:	• UniProtKB - Q3US41 • NCBI - NP_918944.2 • GenelD - 207920

Application Details

Tested Applications:	WB
Application Note:	This protein-A purified antibody has been tested for use western blotting. Specific conditions for reactivity should be optimized by the end user. Expect a band approximately 75.5 kDa in size corresponding to Esrp-1 by western blotting in the appropriate cell lysate or extract.
Assay Dilutions:	All assays should be optimized by the user. Recommended dilutions (if any) may be listed below.
ELISA:	1:750,000
WB:	1:1000

Formulation

Physical State:	Liquid (sterile filtered)
Concentration:	1.0 mg/mL by UV absorbance at 280 nm
Buffer:	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2
Preservative:	0.01% (w/v) Sodium Azide
Stabilizer:	None


Shipping & Handling

Shipping Condition:	Dry Ice
----------------------------	---------

Storage Condition: Store vial at -20° C prior to opening. Aliquot contents and freeze at -20° C or below for extended storage. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use.


Expiration: Expiration date is one (1) year from date of receipt.

Images



Western Blot

Protein expression levels of the splicing factors bound to the putative motifs in pressure-overload cardiac hypertrophy. (A) Heart weight/body weight ratio (HW/BW) at 1 week after Sham or TAC operation (left panel). The relative expression levels of the hypertrophic markers in the mouse models were determined by qRT-PCR (right panel). The differences are statistically significant at * $p < 0.05$, ** $p < 0.01$ and *** $p < 0.001$ (The student t-test). Bars represent means $\pm$ SD ($n = 3$). (B) The relative protein expression levels of the putative splicing factors bound to the identified intronic motifs were measured by Western blotting. Bars represent means $\pm$ SD ($n = 3$). The differences are statistically significant at * $p < 0.05$, ** $p < 0.01$ (The student t-test). Fig 3. PMID: 24552714



Western Blot

Anti-ESRP-1 by western blot shows detection of ESRP-1 seen in lane 2. Lane 1: GFP-transfected-293T cell extracts, Lane 2: ESRP-1-transfected-293T cell extracts, Lane 3: ESRP2-transfected 293T cell lysates. Briefly, each lane contains approximately 5 μ g of lysate. Primary antibody was used at a 1:1000 dilution in PBS-T plus milk, and reacted for 1hr at room temperature. The membrane was washed and reacted with a 1:10,000 dilution of an anti-mouse ECL antibody for 1hr at room temperature. Molecular weight estimation was made by comparison to prestained MW markers.

References

- Peart NJ et al. The global Protein-RNA interaction map of ESRP1 defines a post-transcriptional program that is essential for epithelial cell function. *iScience*. (2022)
- Kim T et al. Pressure-overload cardiac hypertrophy is associated with distinct alternative splicing due to altered expression of splicing factors. *Mol Cells*. (2014)

Disclaimer

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 5199, Limerick, Pennsylvania, USA.