

Datasheet for 500-301-980**H5N1 Antibody VN04-16****Overview**

Description:	Anti-H5 Hemagglutinin of A/Vietnam/1203/04 Influenza Virus (VN04-16) Ascites M (MOUSE) Monoclonal Antibody - 500-301-980
Item No.:	500-301-980
Size:	100 µL
Applications:	ELISA
Reactivity:	Virus
Host Species:	Mouse

Product Details

Background:	Antibody raised against the hemagglutinin (HA) surface glycoprotein of the A/Vietnam/1203/04 (H5N1) influenza virus. Generally referred to as "bird flu", the H5N1 influenza A virus has been documented in poultry and humans across ten Eurasian countries, from Japan in the north to Indonesia in the south. Without immunity, humans would have no protection against H5N1 influenza viruses, which could potentially cause a catastrophic pandemic influenza. This antibody, directed against the HA surface glycoprotein of the A/Vietnam/1203/04 (H5N1) influenza virus, is intended to further our understanding of the mechanisms underlying antigenic variation and evolution of novel variants. The major functions of HA include receptor-binding and fusion activities, but there may also be a structural role for HA in viral particle formation. Following attachment of HA to surface receptors on susceptible cells, the influenza virus enters the cell via endocytosis and membrane fusion.
Synonyms:	mouse anti-H5N1 antibody, mouse anti-Hemagglutinin A antibody, H5HA antibody, Hemagglutinin 5 antibody, H5N1 antibody
Host Species:	Mouse
Clonality:	Monoclonal
Clone ID:	18E1
Format:	ASCITES

Target Details

Gene Name:	HA
-------------------	----

Reactivity:	Virus
Immunogen Type:	Native Protein
Immunogen:	This monoclonal antibody was produced by intraperitoneal immunization of BALB/c mice with concentrated purified virus preparation containing hemagglutinin (HA) protein of influenza A virus [strain A/Vietnam/1203/04 (H5N1)] using the modification of the method described by Kohler and Milstein. Each mouse received two immunizations of 15 µg HA with incomplete Freund's adjuvant, administered 3 week apart.
Purity/Specificity:	This product was clarified from mouse ascitic fluid and is specific for H5 hemagglutinin (HA) protein of influenza A virus [strain A/Vietnam/1203/04 (H5N1)]. VN04-16 monoclonal antibody did not cross-react with influenza viruses of other HA subtypes. This monoclonal antibody reacted with H5N1 influenza virus representatives of different clades and subclades of the H5 HA subtype.
Relevant Links:	• NCBI - 159144921 • UniProtKB - A8UDQ2 • SDS

Application Details

Suggested Applications:	ELISA (Based on references)
Application Note:	This monoclonal antibody can be used for hemagglutination inhibition (HI) assays to provide antigenic characterization of the influenza A viruses of the H5 HA subtype. This monoclonal antibody is suitable for virus neutralization assays (in cell culture and in embryonated chicken eggs), ELISA, immunoprecipitation, immunohistochemistry and western blotting.
Assay Dilutions:	All assays should be optimized by the user. Recommended dilutions (if any) may be listed below.
ELISA:	1:5,000
IHC:	User Optimized
IP:	User Optimized
Neutralization:	User Optimized
WB:	User Optimized

Formulation

Physical State:	Liquid (sterile filtered)
Buffer:	None
Preservative:	0.01% (w/v) Sodium Azide

Stabilizer: None


Shipping & Handling

Shipping Condition: Dry Ice

Storage Condition: Store vial at -20° C prior to opening. Aliquot contents and freeze at -20° C or below for extended storage. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use.

Expiration: Expiration date is one (1) year from date of receipt.

Images




ELISA

Relative reactivity of monoclonal antibody against various VN04 constructs. ELISA plates were coated with 4 µg/ml of each protein in duplicates overnight, blocked and incubated with a 1:5000 dilution of each of the monoclonal antibodies (p/n 500-301-975 through 500-301-980) for 2 h at room temperature (one antibody per plate), followed by a 30 min incubation of a 1:10,000 dilution of HRP-goat anti-mouse IgG for 30 min and developed with TMB. Mean absorbance with standard deviation for replicate wells are reported. Fig 3. PMID: 19654064

Diagram

Schematic representation of the antigenic sites and the epitopes on the globular head of the HA H5 HA molecule. Images were created with RasMol 2.6, and the HA structure was obtained from the Protein Data Bank (PDB accession number 1JSM). Amino acid positions are designated in H3 numbering. Image provided courtesy of Elena Govorkova Ph D.



References

- Song L et al. Superior efficacy of a recombinant flagellin:H5N1 HA globular head vaccine is determined by the placement of the globular head within flagellin. *Vaccine*. (2009)

Disclaimer

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 5199, Limerick, Pennsylvania, USA.