

Datasheet for 600-101-197**APOLIPOPROTEIN E Antibody****Overview**

Description:	Anti-Apolipoprotein E (GOAT) Antibody - 600-101-197
Item No.:	600-101-197
Size:	500 µg
Applications:	WB, Other
Reactivity:	Human, Rat
Host Species:	Goat

Product Details

Background:	Anti Apolipoprotein E Antibody recognizes the gene product of APOE. Apolipoprotein E is found in the chylomicron and IDLs that binds to a specific receptor on liver cells and peripheral cells. Apolipoprotein E is essential for the normal catabolism of triglyceride-rich lipoprotein constituents. This antibody is suitable for cardiovascular research.
Synonyms:	goat anti-Apolipoprotein E3 antibody, Alzheimer disease 2 antibody, Apo-E antibody, Apo E antibody, ApoE antibody, APOEA antibody, Apolipoprotein-E antibody, Apoprotein antibody
Host Species:	Goat
Clonality:	Polyclonal
Format:	IgG

Target Details

Gene Name:	APOE
Reactivity:	Human, Rat
Immunogen Type:	Native Protein
Immunogen:	Apolipoprotein Type E was isolated from human plasma by density gradient centrifugation followed by HPLC purification.

Purity/Specificity: Anti-Apolipoprotein E has been prepared by immunoaffinity chromatography using immobilized antigens followed by extensive cross-adsorption against other apoLipoproteins and human serum proteins to remove any unwanted specificities. Typically less than 1% cross-reactivity against other types of apoLipoprotein was detected by ELISA against purified standards. This antibody reacts with human apoLipoprotein E and has negligible cross-reactivity with Type A-1 A-II, B, C-I, C-II, C-III and J apoLipoproteins. Specific cross-reaction of anti-apoLipoprotein antibodies with antigens from other species has not been determined. Non-specific cross-reaction of anti-apoLipoprotein antibodies with other human serum proteins is negligible.

Relevant Links:

- [600-101-197 SDS](#)
- [UniProtKB - P02649](#)
- [NCBI - NP_000032](#)
- [GenelD - 348](#)

Application Details

Tested Applications:	WB
Suggested Applications:	Other (Based on references)
Application Note:	Apolipoprotein E antibodies have been tested by western blot and used for indirect trapping ELISA for quantitation of antigen in serum using a standard curve, for immunoprecipitation and for western blotting for highly sensitive qualitative analysis.
Assay Dilutions:	All assays should be optimized by the user. Recommended dilutions (if any) may be listed below.
ELISA:	1:4,000 - 1:8,000
IHC:	1:50 - 1:200
IP:	1:100
WB:	1:500 - 1:2,000


Formulation

Physical State:	Liquid (sterile filtered)
Concentration:	1.0 mg/mL by UV absorbance at 280 nm
Buffer:	0.125 M Sodium Borate, 0.075 M Sodium Chloride, 0.005 M EDTA, pH 8.0
Preservative:	0.01% (w/v) Sodium Azide
Stabilizer:	None

Shipping & Handling

Shipping Condition:	Wet Ice
Storage Condition:	Store vial at 4° C prior to opening. This product is stable at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, mix with an equal volume of glycerol, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing.
Expiration:	Expiration date is one (1) year from date of receipt.

Images



Western Blot


Analysis of plasma TC, TG, LDL-C, and HDL-C in the LDLR-KO rabbits.

A: Analysis of plasma apolipoproteins by Western blot. Plasma samples (0.5 μ L) were fractionated on 4–15% SDS-PAGE and transferred to a PVDF membrane probed with Abs anti ApoA-I, ApoB, and ApoE as described in the Materials and Methods section.

B: Agarose gel electrophoresis of plasma lipoproteins. 4 μ L of plasma was loaded in each well, fractionated on 1% agarose gel, and stained with fat red 7 B for neutral lipids. Lipoprotein migration positions are indicated by arrows.

C: Analysis of plasma lipoprotein profiles by fast-protein liquid chromatography (FPLC). Black, red, green, and blue lines showed the fractions of WT, L9, L12, and L17 respectively; VLDL, very low-density lipoprotein; LDL, low-density lipoprotein; HDL, high-density lipoprotein.

Fig 3. PMID: 30243490




Western Blot

CBS-KO rabbits are hyperlipidemic. a Plasma TG, TC, HDL-C, LDL-C and b Hcy levels in 6 weeks-old WT rabbits (n = 6). Mean $\pm$ SD, *P < 0.05. c Immunoblot showing levels of plasma apolipoproteins in the indicated groups. WT-NC: wild type rabbits fed normal chow. WT-SB: wild type rabbits fed normal chow supplemented with the vitamin B and betaine complex. C2, C8: CBS-KO rabbits on a normal chow diet. C9, C14, C15, C16: CBS-KO rabbits fed normal chow supplemented with the vitamin B and betaine complex. Fig 3. PMID: 33054837

Western Blot

Western Blot of goat Anti Apolipoprotein E Antibody. Lane 1: rat glial whole cell lysate. Load: 40 μ g per lane. Primary antibody: Apolipoprotein E antibody at 1:1000 for overnight at 4°C. Secondary antibody: HRP goat secondary antibody at 1:200 for 45 min at RT. Block: 3%BLOTTO overnight at 4°C. Predicted/Observed size: 36 kDa, 35-36 kDa for APO-E. MW markers arrows at 25kDa and 37kDa. Other band(s): ~55kDa.



References

- Zhang T et al. Hyperhomocysteinemia and Dyslipidemia in point mutation G307S of cystathionine β -synthase-deficient rabbit generated using CRISPR/Cas9. *Lipids Health Dis.* (2020)
- Lu et al. Spontaneous severe hypercholesterolemia and atherosclerosis lesions in rabbits with deficiency of low-density lipoprotein receptor (LDLR) on exon 7. *EBioMedicine* (2018)
- Berard AR et al. Quantification of the host response proteome after herpes simplex virus type 1 infection. *J Proteome Res.* (2015)
- Berard AR et al. Comparative proteomic analyses of two reovirus T3D subtypes and comparison to T1L identifies multiple novel proteins in key cellular pathogenic pathways. *Proteomics.* (2015)

Disclaimer

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 5199, Limerick, Pennsylvania, USA.