

Datasheet for 600-401-445**IRS1 phospho S307 Antibody****Overview**

Description:	Anti-IRS1 pS307 (RABBIT) Antibody - 600-401-445
Item No.:	600-401-445
Size:	100 µg
Applications:	ELISA, WB
Reactivity:	Human
Host Species:	Rabbit

Product Details

Background:	Insulin Receptor Substrate 1 (IRS1) acts as a signaling molecule for IL-4, insulin and insulin-like growth factor-I (IGF-I) receptors. When phosphorylated by the insulin receptor, IRS1 binds specifically to various cellular proteins containing SH2 domains such as phosphatidylinositol 3-kinase p85 subunit or GRB2. When bound, IRS1 typically activates phosphatidylinositol 3-kinase p85 subunit. IRS1 interacts with both the NPXY motif of tyrosine-phosphorylated IGF1R and the INSR through the PTB domain. Serine phosphorylation of IRS1 is a mechanism for insulin resistance. Diseases associated with IRS1 include Diabetes Mellitus, Noninsulin-Dependent and Mixed Cell Adenoma. Anti-IRS1pS307 is useful for researchers interested in protein kinase binding, IL-2 pathways, and RET signaling.
Synonyms:	rabbit anti-IRS1 pS307 antibody, IRS-1, IRS 1, HIRS 1 antibody, HIRS-1 antibody, Insulin Receptor Substrate 1 antibody
Host Species:	Rabbit
Clonality:	Polyclonal
Format:	IgG

Target Details

Gene Name:	IRS1
Reactivity:	Human
PTM Specificity:	Phosphorylation
Immunogen Type:	Conjugated Peptide

Immunogen:	This affinity purified antibody was prepared from whole rabbit serum produced by repeated immunizations with a synthetic phospho-peptide corresponding to an internal region near amino acids 290-320 of human IRS1 protein.
Purity/Specificity:	This affinity-purified antibody is directed against the phosphorylated form of human IRS1 protein at the pS307 residue. The product was affinity purified from monospecific antiserum by immunoaffinity purification. Antiserum was first purified against the phosphorylated form of the immunizing peptide. The resultant affinity purified antibody was then cross-adsorbed against the non-phosphorylated form of the immunizing peptide. Reactivity occurs against human IRS1 pS307 protein and the antibody is specific for the phosphorylated form of the protein. Reactivity with non-phosphorylated human IRS1 is minimal by ELISA. A BLAST analysis was used to suggest cross reactivity with IRS1 from human, mouse, rat, dog and vervet monkey based on 100% homology with the immunizing sequence. Partial reactivity is also expected against IRS1 from pig (94%), bovine (94%) and chicken (88%) sources. Reactivity of this antibody with IRS1 from other species is unknown.
Relevant Links:	• UniProtKB - P35568 • NCBI - NP_005535.1 • GenelD - 3667

Application Details

Tested Applications:	ELISA, WB
Application Note:	This affinity purified antibody has been tested for use in ELISA and western blot. Specific conditions for reactivity should be optimized by the end user. Expect a band approximately 130 kDa in size corresponding to phosphorylated IRS1 protein by western blotting in the appropriate cell lysate or extract. Minimal reactivity is observed against the non-phosphorylated form of the immunizing peptide. This antibody is phospho specific for pS307 of IRS1 protein.
Assay Dilutions:	All assays should be optimized by the user. Recommended dilutions (if any) may be listed below.
ELISA:	1:2,000 - 1:10,000
IHC:	User Optimized
WB:	1:250 - 1:1,500

Formulation

Physical State:	Liquid (sterile filtered)
Concentration:	1.0 mg/mL by UV absorbance at 280 nm
Buffer:	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2

Preservative: 0.01% (w/v) Sodium Azide

Stabilizer: None


Shipping & Handling

Shipping Condition: Dry Ice

Storage Condition: Store vial at -20° C prior to opening. Aliquot contents and freeze at -20° C or below for extended storage. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use.

Expiration: Expiration date is one (1) year from date of receipt.

Images



Western Blot

Western blot using Rockland's Affinity Purified anti-IRS1 pS307 antibody shows detection of a band at ~180 kDa believed to represent phosphorylated IRS1 (arrowhead). Lane 1 shows staining of human 293 cell lysate (p/n W09-000-365). Lane 2 shows staining of 293 cell lysate prepared from cells serum-starved for 18 h followed by treatment with 5 µg/ml of anisomycin for 30 min. The pronounced staining of the band at 180 kDa is not seen when the antibody was pre-incubated with immunizing peptide prior to reaction (data not shown). The identity of the intensely reactive bands at ~70 kD in both lane 1 and 2 is unknown, although these bands were also competed out by pre-incubation with the immunizing peptide. Approximately 25 µg of each lysate was separated on a 4-20% Tris-Glycine gel by SDS-PAGE and transferred onto nitrocellulose. After blocking with 5% goat serum, 0.5% BLOTTO in PBS, the membrane was probed with the primary antibody diluted to 1:250. Reaction occurred overnight at 4° C followed by washes and reaction with a 1:10,000 dilution of IRDye™800 conjugated Gt-a-Rabbit IgG [H&L] MX (611-132-122) for 45 min at room temperature (800 nm channel, green). Molecular weight estimation was made by comparison to prestained MW markers in lane M (700 nm channel, red). IRDye™800 fluorescence image was captured using the Odyssey® Infrared Imaging System developed by LI-COR. IRDye is a trademark of LI-COR, Inc. Other detection systems will yield similar results.

Disclaimer

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 5199, Limerick, Pennsylvania, USA.