

Datasheet for 600-401-448**RFA2 Antibody****Overview**

Description:	Anti-RFA2 pan reactive (RABBIT) Antibody - 600-401-448
Item No.:	600-401-448
Size:	100 µg
Applications:	ELISA, WB
Reactivity:	Yeast
Host Species:	Rabbit

Product Details

Background:	RFA2 binds to single-stranded sequences participating in DNA replication in addition to those mediating transcriptional repression (URS1) and activation (CAR1). It stimulates the activity of a cognate strand exchange protein (SEP1). It cooperates with T-AG and DNA topoisomerase I to unwind template DNA containing the simian virus 40 origin of DNA replication.
Synonyms:	rabbit anti-RFA2 antibody, RFA-2, RFA 2, Replication factor A protein 2, RF-A protein 2, Replication protein A 36 kDa subunit, DNA-binding protein BUF1, BUF-1, YNL312W, N0368
Host Species:	Rabbit
Clonality:	Polyclonal
Format:	IgG

Target Details

Gene Name:	RFA2
Reactivity:	Yeast
Immunogen Type:	Conjugated Peptide
Immunogen:	This affinity purified antibody was prepared from whole rabbit serum produced by repeated immunizations with a synthetic peptide corresponding to an internal region near aa 110-135 (of 273) of yeast DNA Replication Factor A (RFA2). The RFA complex consists of a heterotrimer of RFA1, RFA2 and RFA3, 69, 36 and 13 kda chains respectively. Yeast RFA2 is phosphorylated at the G₁/S transition and is dephosphorylated at mitosis.

Purity/Specificity: This affinity purified antibody is directed against yeast RFA2. The product was affinity purified from monospecific antiserum by immunoaffinity purification. Antiserum was purified against the phosphorylated form of the immunizing peptide and then cross adsorbed against the non-phosphorylated form of the peptide. This antibody consists of antibodies eluted from the non-phospho peptide column.

Relevant Links:

- [UniProtKB - P26754](#)
- [NCBI - 132475](#)
- [GeneID - 855404](#)

Application Details

Tested Applications: ELISA, WB

Application Note: This affinity-purified antibody was tested for use in western blotting against RFA2 containing cell lysates and by ELISA against the immunizing peptide. Reactivity in other immunoassays is unknown.

Assay Dilutions: All assays should be optimized by the user. Recommended dilutions (if any) may be listed below.

ELISA: 1:1000 - 1:3000

IP: 1:100

WB: 1:100 - 1:500

Formulation

Physical State: Liquid (sterile filtered)

Concentration: 0.12 mg/mL by UV absorbance at 280 nm

Buffer: 0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2

Preservative: 0.01% (w/v) Sodium Azide

Stabilizer: None


Shipping & Handling

Shipping Condition: Dry Ice

Storage Condition: Store vial at -20° C prior to opening. Aliquot contents and freeze at -20° C or below for extended storage. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use.

Expiration: Expiration date is one (1) year from date of receipt.

Images



Western Blot

This product was assayed by western blotting against RFA2 containing cell lysates. Blot incubated for 1 hr at RT with affinity purified Rabbit anti RFA2 pan reactive Lot 13229 at 1.0 ug/ml followed by 1:5000 dilution of IRDye800 Sheep anti Rabbit IgG (lot 13218) for 1 hr at RT.

Disclaimer

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 5199, Limerick, Pennsylvania, USA.