

Datasheet for 600-401-466**KLH Antibody****Overview**

Description:	Anti-KLH (RABBIT) Antibody - 600-401-466
Item No.:	600-401-466
Size:	1 mg
Applications:	ELISA, WB
Reactivity:	KLH
Host Species:	Rabbit

Product Details

Background:	KLH Antibody detects Keyhole Limpet Hemocyanin protein which is a very large, heterogeneous glycosylated protein consisting of subunits with a molecular weight of 350,000 and 390,000. KLH readily forms aggregates with molecular weights of 4,500,000-13,000,000 daltons. Each domain of a KLH subunit contains two copper atoms that together bind a single oxygen molecule (O ₂). When oxygen is bound to hemocyanin, the molecule takes on a distinctive transparent, opalescent blue color. The KLH protein is potentially immunogenic yet safe in humans and is therefore highly prized as a vaccine carrier protein. The large and highly glycosylated KLH protein cannot be reproduced synthetically. It is available only as a purified biological product from the Keyhole Limpet <i>Megathura crenulata</i> . KLH is being tested in a variety of cancer vaccines, including non-Hodgkins lymphoma, cutaneous melanoma, breast and bladder cancer.
Synonyms:	rabbit Anti-KLH antibody, Keyhole Limpet Hemocyanin
Host Species:	Rabbit
Clonality:	Polyclonal
Format:	IgG

Target Details

Reactivity:	KLH
Immunogen Type:	Native Protein
Immunogen:	KLH (Keyhole Limpet Hemocyanin)

Purity/Specificity:	Anti-KLH antibody was prepared from monospecific antiserum by immunoaffinity chromatography using KLH coupled to agarose beads followed by solid phase adsorption(s) to remove any unwanted reactivities. Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-Rabbit Serum and KLH.
----------------------------	---

Application Details

Tested Applications:	ELISA, WB
Application Note:	Anti-KLH Antibody has been tested for western blotting to detect a single band of the expected apparent molecular weight and for ELISA. Researchers should determine optimal titers for applications that are not stated below.
Assay Dilutions:	All assays should be optimized by the user. Recommended dilutions (if any) may be listed below.
ELISA:	1:10,000 - 1:50,000
IHC:	1:500 - 1:2,500
WB:	1:1,000 - 1:5,000


Formulation

Physical State:	Liquid (sterile filtered)
Concentration:	1.0 mg/mL by UV absorbance at 280 nm
Buffer:	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2
Preservative:	0.01% (w/v) Sodium Azide
Stabilizer:	None

Shipping & Handling

Shipping Condition:	Wet Ice
Storage Condition:	Store KLH antibody at 4° C prior to opening. Aliquot contents and freeze at -20° C or below for extended storage. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use.
Expiration:	Expiration date is one (1) year from date of receipt.

Images

**Western Blot**

Western Blot of Rabbit anti-KLH antibody. Lane 1: Molecular Weight Marker. Lane 2: KLH. Load: 10 µg. Primary antibody: Anti-KLH at 1:1000 in 1% BSA-TBS, overnight at 4°C. Secondary antibody: HRP rabbit secondary antibody at 1:40,000 for 45 min at RT. Block: 1% BSA-TTBS for 30 min RT. Predicted/Observed size: ~350-390 kDa for KLH.

Disclaimer

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 5199, Limerick, Pennsylvania, USA.