

Datasheet for 600-401-485**Fbp5B Antibody****Overview**

Description:	Anti-F-Box Protein 43 (Fbp5B) (RABBIT) Antibody - 600-401-485
Item No.:	600-401-485
Size:	100 µg
Applications:	ELISA, WB
Reactivity:	Human, Mouse, Rat, Chimpanzee
Host Species:	Rabbit

Product Details

Background:	Fbp5B (also called F-Box protein 43 or Fbx043) is a relatively new member of the F-box protein family that is characterized by an approximately 40 amino acid motif, the F-box. The F-box proteins constitute one of the four subunits of the ubiquitin protein ligase complex called SCFs (SKP1-cullin-F-box), which function in phosphorylation-dependent ubiquitination. The F-box proteins are divided into 3 classes: Fbws containing WD-40 domains, Fbls containing leucine-rich repeats, and Fbxs containing either different protein-protein interaction modules or no recognizable motifs. The protein encoded by this gene belongs to the Fbxs class.
Synonyms:	rabbit anti-F-Box Protein 43/Fbp5B antibody, rabbit anti-F-Box Protein 43 antibody, rabbit anti-Fbp5B antibody, F-box only protein 43, Endogenous meiotic inhibitor 2, FBXO43, EMI2
Host Species:	Rabbit
Clonality:	Polyclonal
Format:	IgG

Target Details

Gene Name:	FBXO43
Reactivity:	Human, Mouse, Rat, Chimpanzee
Immunogen Type:	Conjugated Peptide

Immunogen:	This affinity purified antibody was prepared from whole rabbit serum produced by repeated immunizations with a synthetic peptide corresponding to a region approximately 60 residues downstream of the amino terminal end of human Fbp5B protein. Human Fbp5B-2, 466 aa, predicted MW 51.3 kDa
Purity/Specificity:	This affinity-purified antibody is directed against human Fbp5B protein. The product was affinity purified from monospecific antiserum by immunoaffinity purification. A BLAST analysis was used to suggest cross reactivity with Fbp5B protein from human and chimpanzee based on 100% homology with the immunizing sequence. Reactivity against homologues from other sources is not known.
Relevant Links:	• NCBI - 117606351 • UniProtKB - Q4G163 • NCBI - Q4G163.3 • GeneID - 286151

Application Details

Tested Applications:	ELISA, WB
Application Note:	This affinity-purified antibody has been tested for use in ELISA and western blot. Specific conditions for reactivity should be optimized by the end user. Expect a band approximately 105 kDa in size corresponding to Fbp5B protein by western blotting in the appropriate cell lysate or extract. Predicted molecular weight is 51 kDa.
Assay Dilutions:	All assays should be optimized by the user. Recommended dilutions (if any) may be listed below.
ELISA:	1:10,000 - 1:50,000
IP:	1:100
WB:	1:500 - 1:3,000


Formulation

Physical State:	Liquid (sterile filtered)
Concentration:	2.21 mg/mL by UV absorbance at 280 nm
Buffer:	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2
Preservative:	0.01% (w/v) Sodium Azide
Stabilizer:	None

Shipping & Handling

Shipping Condition:	Dry Ice
Storage Condition:	Store vial at -20° C prior to opening. Aliquot contents and freeze at -20° C or below for extended storage. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use.
Expiration:	Expiration date is one (1) year from date of receipt.

Images



Western Blot

Western blot using Rockland's affinity purified anti-Fbp5b antibody. Panels A and B show two identical membranes containing marker, in vitro translated Fbp5b (Fbp5b ivt), and a whole cell extract from 293T cells loaded in duplicate). Panel A shows antibody reactivity. Panel B shows antibody reactivity after first pre-incubating the antibody with the immunizing peptide. This pre-incubation greatly diminishes specific band recognition. Peptide inhibition completely removes the in vitro translated band and greatly reduces the endogenous band from 293T cells. When this experiment was repeated using a mock in vitro translation, no band was detected (data not shown). Personal Communication, Daniele Guardavaccaro, NYU Cancer Institute

Disclaimer

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 5199, Limerick, Pennsylvania, USA.