

Datasheet for 600-401-487**APG7 Antibody****Overview**

Description:	Anti-APG7 (RABBIT) Antibody - 600-401-487
Item No.:	600-401-487
Size:	100 µg
Applications:	ELISA, WB
Reactivity:	Human
Host Species:	Rabbit

Product Details

Background:	Apg7 is the human gene product similar in homology to Pichia pastoris GSA7 and Saccharomyces cerevisiae Apg7. In the yeast, the protein appears to be required for fusion of peroxisomal and vacuolar membranes. The protein shows homology to the ATP-binding and catalytic sites of the E1 ubiquitin activating enzymes.
Synonyms:	rabbit anti-APG7 Antibody, APG-7, APG 7, APG7 autophagy 7-like antibody, APG7 like antibody, APG7L antibody, ATG 7 antibody, ATG-7 antibody, ATG7 autophagy related 7 homolog antibody
Host Species:	Rabbit
Clonality:	Polyclonal
Format:	IgG

Target Details

Gene Name:	ATG7
Reactivity:	Human
Immunogen Type:	Conjugated Peptide
Immunogen:	This affinity purified antibody was prepared from whole rabbit serum produced by repeated immunizations with a synthetic peptide corresponding to an internal region near aa 550-575 of Human Apg7 protein.

Purity/Specificity: This is an affinity purified antibody produced by immunoaffinity chromatography using the immunizing peptide after immobilization to a solid phase. Reactivity occurs against human Apg7 protein. However, BLAST analysis indicates 100% homology (12/12) for this protein from human, chimpanzee, mouse, rat, dog, and chicken sources. The core amino acids of the immunogen, G-D-S-T-R-D-R-T are 100% identical (8/8) in D. melanogaster. Reactivity with Apg7 proteins from other sources is not known.

Relevant Links:

- [UniProtKB - O95352](#)
- [NCBI - NP_001129503.2](#)
- [GeneID - 10533](#)

Application Details

Tested Applications: ELISA, WB

Application Note: This affinity purified antibody has been tested for use in ELISA and by western blot. Specific conditions for reactivity should be optimized by the end user. Expect a band approximately 78 kDa in size corresponding to human Apg7 by western blotting in the appropriate cell lysate or extract. Human Apg7 has been reported to form a homodimer at approximately 160 kDa.

Assay Dilutions: All assays should be optimized by the user. Recommended dilutions (if any) may be listed below.

ELISA: 1:300,000 - 1:400,000

IF: User Optimized

WB: 1:500 - 1:2,000

Formulation

Physical State: Liquid (sterile filtered)

Concentration: 1.0 mg/mL by UV absorbance at 280 nm

Buffer: 0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2

Preservative: 0.01% (w/v) Sodium Azide

Stabilizer: None


Shipping & Handling

Shipping Condition: Dry Ice

Storage Condition: Store vial at -20° C prior to opening. Aliquot contents and freeze at -20° C or below for extended storage. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use.


Expiration: Expiration date is one (1) year from date of receipt.

Images




Western Blot

(C) HeLa cells treated with or without 100mM trehalose for 24h were analyzed for Atg7 levels by immunoblotting with anti-Atg7 (p/n 600-401-487) antibody (i) and densitometry analysis relative to Tubulin (ii). $p=0.4315$. Fig. S5. PMID: 17182613




Western Blot

The effect of SMERs10, 18 and 28 on mTOR activity, Beclin-1/Atg6, Atg5, Atg7, Atg12, Atg5-Atg12 conjugation and proteasome activity. (f) HeLa cells treated with DMSO (control) or with 47 μM SMER10, 43 μM SMER18 or 47 μM SMER28 for 24 h, were analyzed for Atg7 (f) levels by immunoblotting with anti-Atg7 (p/n 600-401-487) antibodies. Supplementary Figure 3. PMID: 17486044




Western Blot

D, BMDM from *Atg7^{f/f}* or *Atg7^Δ* mice were left untreated, primed with LPS, or primed with LPS and then treated with ATP (5 mM) for 30 min (left panel) or 0–60 min (right panel). Fig 1. PMID: 24178293



Western Blot

ALIS formation and clearance are independent of ATG7. BMMs from wild type (*ATG7^{f/f}*) or *ATG7* knock-out (*ATG7^{f/f}Cre⁺*) mice were infected with BCG for 1 h followed by culture in regular medium for 7 h (BCG8h) or 35 h (BCG36h). In some experiments, BMMs were infected with BCG for 1 h followed by culture in regular medium for 7 h and then treated for 4 h with 25 μM cycloheximide either with or without 200 nM of bafilomycin A. Cell lysates were subjected to immunoblot analysis using antibodies against ATG7, p62, LC3B, or β-actin (A). Fig 5. PMID: 22518844



Western Blot

(D) Whole cell lysates of HeLa-R30 cells were collected 48 h after transfection of siRNA targeted against ATG7 (siATG7) or non-targeted siControl and immunoblotted for ATG7, LC3 or β-tubulin. Greater than 80% knockdown was achieved with siATG7 when compared with siControl. LC3-I accumulation was noted in cells with ATG7 knockdown. Fig 4. PMID: 22954701

References

- Fattah et al. Critical role for IL-18 in spontaneous lung inflammation caused by autophagy deficiency. *Journal of Immunology* (2015)
- Bhattacharya A et al. Deficiency of autophagy in dendritic cells protects against experimental autoimmune encephalomyelitis. *J Biol Chem.* (2014)
- Shin JN et al. Inflammasome activation by altered proteostasis. *J Biol Chem.* (2013)
- Fung C et al. EGFR tyrosine kinase inhibition induces autophagy in cancer cells. *Cancer Biol Ther.* (2012)
- Liu XD et al. Transient aggregation of ubiquitinated proteins is a cytosolic unfolded protein response to inflammation and endoplasmic reticulum stress. *J Biol Chem.* (2012)
- Sarkar S et al. Small molecules enhance autophagy and reduce toxicity in Huntington's disease models. *Nat Chem Biol.* (2007)
- Sarkar S et al. Trehalose, a novel mTOR-independent autophagy enhancer, accelerates the clearance of mutant huntingtin and α -synuclein. *J Biol Chem.* (2007)

Disclaimer

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 5199, Limerick, Pennsylvania, USA.