

Datasheet for 600-401-A89S**Pogz Antibody****Overview**

Description:	Anti-Pogz (RABBIT) Antibody - 600-401-A89S
Item No.:	600-401-A89S
Size:	25 µL
Applications:	ELISA, WB, ChIP, IP
Reactivity:	Mouse
Host Species:	Rabbit

Product Details

Background:	This antibody is designed, produced, and validated as part of a collaboration between Rockland and the National Cancer Institute (NCI) and is suitable for Cancer, Immunology and Nuclear Signaling research. Pogz is a zinc-finger protein containing at least 8 C2H2 zinc fingers, a CENP-B (centromere protein-B) domain, and a DDE domain found in the bacterial transposase/retroviral integrase family. Members of this family of integrases share structural homologies and show similar catalytic activity to the RAG1 and RAG2 proteins that initiate V(D)J recombination in lymphoid progenitors. Mouse Pogz has a human homolog with a base pair similarity of 90% and amino acid similarity of 93% for the full length protein. In humans at least 3 variants have been predicted and they are all similar at the N-terminus. A number of expressed sequence tags (ESTs) cloned from murine undifferentiated ES cells and Lin-/c-Kit+/Sca-1+ hematopoietic stem cell cDNA libraries correspond to the Pogz gene, further underlining an important function for this gene in stem cells. In mice, deletion of the gene in all tissues early in embryogenesis has been shown to be lethal.
Synonyms:	rabbit anti-POGZ antibody, pogo transposable element with ZNF domain antibody, Suppressor of hairy wing homolog 5, Zinc finger protein 280E, Zinc finger protein 635, SUHW5, ZNF280E, ZNF635
Host Species:	Rabbit
Clonality:	Polyclonal
Format:	IgG

Target Details

Gene Name:	Pogz
-------------------	------

Reactivity:	Mouse
Immunogen Type:	Conjugated Peptide
Immunogen:	This affinity purified antibody was prepared from whole rabbit serum produced by repeated immunizations with a synthetic peptide corresponding to a region near the N-terminus of mouse Pogz protein.
Purity/Specificity:	This product was affinity purified from monospecific antiserum by immunoaffinity chromatography. This antibody is known to react with mouse Pogz protein. A BLAST analysis was used to suggest cross-reactivity with Pogz from human, dog, short-tailed opossum, cattle, rat, chimpanzee, macaque, olive baboon, and chicken sources based on 100% homology with the immunizing sequence. Reactivity with Pogz from other sources has not been determined.
Relevant Links:	• UniProtKB - Q8BZH4 • NCBI - Q8BZH4.2 • GeneID - 229584

Application Details

Tested Applications:	ELISA, WB
Suggested Applications:	ChIP, IP (Based on references)
Application Note:	This affinity purified antibody has been tested for use in ELISA and western blotting. Specific conditions for reactivity should be optimized by the end user. Expect a band approximately 155 kDa in size corresponding to Pogz by western blotting in the appropriate cell lysate or extract.
Assay Dilutions:	All assays should be optimized by the user. Recommended dilutions (if any) may be listed below.
ELISA:	1:150,000
WB:	1:3,000


Formulation

Physical State:	Liquid (sterile filtered)
Concentration:	1.6 mg/mL by UV absorbance at 280 nm
Buffer:	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2
Preservative:	0.01% (w/v) Sodium Azide
Stabilizer:	None

Shipping & Handling


Shipping Condition:	Dry Ice
Storage Condition:	Store vial at -20° C or below prior to opening. This vial contains a relatively low volume of reagent (25 µL). To minimize loss of volume dilute 1:10 by adding 225 µL of the buffer stated above directly to the vial. Recap, mix thoroughly and briefly centrifuge to collect the volume at the bottom of the vial. Use this intermediate dilution when calculating final dilutions as recommended below. Store the vial at -20°C or below after dilution. Avoid cycles of freezing and thawing.
Expiration:	Expiration date is three (3) months from date of receipt.

Images




ChIP

(A) Heatmap of ChIP-seq data showing the enrichment of POGZ, SETDB1, and TRIM28 on the indicated ERVs in control and *Pogz*^{-/-} ESCs.
(B) Genomic view of ChIP-seq data for POGZ, SETDB1, and TRIM28 on the Dux locus. Fig 6. PMID: 37494184



Immunoprecipitation

(C) IP images showing that endogenous POGZ pulled down TRIM28 in ESCs (n = 3 repeats).
(D) IP images showing that POGZ pulled down SETDB1 in ESCs (n = 3).
(E) coIP images showing that FLAG-POGZ interacts with MYC-TRIM28 (n = 3).
(F) coIP images showing that FLAG-POGZ interacts with MYC-SETDB1 (n = 3). Fig 6. PMID: 37494184



Western Blot

Western blot using Rockland's affinity purified anti-Pogz antibody shows detection of Pogz protein ~155kDa (arrowhead) in adult mouse thymus and spleen tissue lysate (Panel A). The lower molecular weight bands may be cross reactive proteins. Pre-incubation of antibody with the immunizing peptide blocks specific antibody reactivity (Panel B). Primary antibody was used at 1:20,000. Personal Communication, K.O. Gudmundsson and J. Keller, NCI, Frederick, MD.

References

- Sun X et al. POGZ suppresses 2C transcriptional program and retrotransposable elements. *Cell Rep.* (2023)

Disclaimer

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 5199, Limerick, Pennsylvania, USA.