

Datasheet for 600-401-AU9**DENR Antibody****Overview**

Description:	Anti-DENR (RABBIT) Antibody - 600-401-AU9
Item No.:	600-401-AU9
Size:	100 µg
Applications:	ELISA, IF, IHC, WB
Reactivity:	Human, Mouse, Rat
Host Species:	Rabbit

Product Details

Background:	DENR (density-regulated protein), also designated DRP or smooth muscle cell associated protein 3 (SMAP-3), is a 198 amino acid protein that may be involved in the translation of target mRNAs by scanning and recognition of the initiation codon (1). It plays a role in the modulation of the translational profile of a subset of cancer-related mRNAs when recruited to the translational initiation complex by the oncogene MCTS1 (2,3). DENR contains a SUI1 domain and is found in a variety of tissues with highest levels present in skeletal and cardiac muscle. It is up-regulated in ovarian and breast cancer cells by ERBB2 overexpression (4).
Synonyms:	DENR Antibody, DRP, DRP1, SMAP-3, H14, Density-regulated protein, Protein DRP1, DRP
Host Species:	Rabbit
Clonality:	Polyclonal
Format:	IgG

Target Details

Gene Name:	DENR
Reactivity:	Human, Mouse, Rat
Immunogen Type:	Conjugated Peptide
Immunogen:	Anti-DENR antibody was prepared from whole rabbit serum produced by repeated immunizations with an 18 amino acid peptide near the internal region of human DENR.

Purity/Specificity: Anti-DENR antibody was affinity purified from monospecific antiserum by immunoaffinity chromatography. DENR antibody is human, mouse and rat reactive.

Relevant Links:

- [UniProtKB - O43583](#)
- [GeneID - 8562](#)
- [NCBI - NP_003668](#)

Application Details

Tested Applications: ELISA, IF, IHC, WB

Application Note: Anti-DENR Antibody has been tested for use in ELISA, Western Blotting, Immunohistochemistry and Immunofluorescence. Specific conditions for reactivity should be optimized by the end user. Expect a band at approximately 22 kDa in Western Blots of specific cell lysates and tissues.

Assay Dilutions: All assays should be optimized by the user. Recommended dilutions (if any) may be listed below.

ELISA: 1:10,000 - 1:20,000

IF: 20 µg/mL

IHC: 5 µg/mL

WB: 1-2 µg/mL

Formulation

Physical State: Liquid (sterile filtered)

Concentration: 1 mg/mL by UV absorbance at 280 nm

Buffer: 0.01 M Sodium Phosphate, 0.25 M Sodium Chloride, pH 7.2

Preservative: 0.02% (w/v) Sodium Azide

Stabilizer: None


Shipping & Handling

Shipping Condition: Wet Ice

Storage Condition: Antibody can be stored at 4°C up to one year. Antibodies should not be exposed to prolonged high temperatures.


Expiration: Expiration date is one (1) year from date of receipt.

Images




Immunohistochemistry

Immunohistochemistry of Rabbit anti-DENR antibody.
 Tissue: rat brain tissue. Primary antibody: DENR antibody at 5 µg/mL. Secondary antibody: Peroxidase rabbit secondary antibody at 1:5,000. Staining: DENR as precipitated brown signal.



Immunofluorescence Microscopy

Immunofluorescence Microscopy of Rabbit anti-DENR antibody. Tissue: rat brain. Primary antibody: DENR antibody at 20 µg/mL. Secondary antibody: Fluorescein rabbit secondary antibody at 1:5,000. Staining: DENR as red fluorescent signal.



Western Blot

Western Blot of Rabbit anti-DENR antibody. Lane A: human brain tissue lysate. Primary antibody: DENR antibody at 1 µg/mL overnight at 4°C. Secondary antibody: Goat anti-Rabbit HRP secondary antibody. Block: 5% BLOTTO. Predicted/Observed size: 22 kDa, 22 kDa for DENR.

Disclaimer

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 5199, Limerick, Pennsylvania, USA.