

Datasheet for 600-401-AX1**DPF3 Antibody****Overview**

Description:	Anti-DPF3 (RABBIT) Antibody - 600-401-AX1
Item No.:	600-401-AX1
Size:	100 µg
Applications:	ELISA, IF, IHC, WB
Reactivity:	Human, Mouse
Host Species:	Rabbit

Product Details

Background:	The DPF3 protein, also known as Cerd4, is a member of the d4 gene family of transcription modulators that also includes DPF1/Neud4 and DPF2/Requiem (1). DPF3 has been shown to be a epigenetic key factor for heart and muscle development and can bind to methylated and acetylated lysine residues of histone 3 and 4, suggesting that DPF3 may play a role in recruiting chromatin remodeling complexes to acetylated histones (2). Two isoforms of DPF3, DPF3a and DPF3b, are required as transcriptional co-activators in SWI/SNF complex-dependent activation of the NF-kappaB RelA/p50 heterodimer (3).
Synonyms:	DPF3 Antibody, CERD4, BAF45C, CERD4, Zinc finger protein DPF3, BRG1-associated factor 45C
Host Species:	Rabbit
Clonality:	Polyclonal
Format:	IgG

Target Details

Gene Name:	DPF3
Reactivity:	Human, Mouse
Immunogen Type:	Conjugated Peptide
Immunogen:	Anti-DPF3 antibody was prepared from whole rabbit serum produced by repeated immunizations with an 18 amino acid peptide near the C-terminus of human DPF3

Purity/Specificity:	Anti-DPF3 antibody was affinity purified from monospecific antiserum by immunoaffinity chromatography. DPF3 antibody is human and mouse reactive. Multiple isoforms of DPF3 are known to exist.
Relevant Links:	• UniProtKB - Q92784 • GeneID - 8110 • NCBI - NP_036206

Application Details

Tested Applications:	ELISA, IF, IHC, WB
Application Note:	Anti-DPF3 Antibody has been tested for use in ELISA, Western Blotting, Immunohistochemistry and Immunofluorescence. Specific conditions for reactivity should be optimized by the end user. Expect a band at approximately 43 kDa in Western Blots of specific cell lysates and tissues.
Assay Dilutions:	All assays should be optimized by the user. Recommended dilutions (if any) may be listed below.
ELISA:	1:10,000 - 1:20,000
IF:	20 µg/mL
IHC:	5 µg/mL
WB:	1-2 µg/mL

Formulation


Physical State:	Liquid (sterile filtered)
Concentration:	1 mg/mL by UV absorbance at 280 nm
Buffer:	0.01 M Sodium Phosphate, 0.25 M Sodium Chloride, pH 7.2
Preservative:	0.02% (w/v) Sodium Azide
Stabilizer:	None

Shipping & Handling

Shipping Condition:	Dry Ice
Storage Condition:	Store vial at -20° C prior to opening. Aliquot contents and freeze at -20° C or below for extended storage. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use.

Expiration:


Expiration date is one (1) year from date of receipt.

Images**Western Blot**

Western blot of DPF3.

Load: mouse brain tissue lysate.


Primary Antibody: DPF3 antibody at (A) 1 and (B) 2 μ g/ml.

**Immunohistochemistry**

Immunohistochemistry of DPF3.

Tissue: human brain tissue.

Primary Antibody: DPF3 antibody at 5 μ g/mL.

**Immunofluorescence Microscopy**

Immunofluorescence of DPF3.

Tissue: human brain tissue.

Primary Antibody: DPF3 antibody at 20 µg/mL.

Disclaimer

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 5199, Limerick, Pennsylvania, USA.