

Datasheet for 600-401-CH4**LAMTOR2 Antibody****Overview**

Description:	Anti-LAMTOR2 (RABBIT) Antibody - 600-401-CH4
Item No.:	600-401-CH4
Size:	100 µg
Applications:	ELISA, IF, IHC, WB
Reactivity:	Human, Mouse, Rat
Host Species:	Rabbit

Product Details

Background:	The late endosomal/lysosomal adaptor MAPK and MTOR activator 2 (LAMTOR2) protein belongs to the LAMTOR family of proteins, and together with LAMTOR3 and the MAPK1 and ERK kinase 1 (MEK1) localizes to late endosomes where it is required for the efficient activation of ERK signaling (1,2). This complex is involved in the regulation of late endosomal traffic and cellular proliferation (3) and plays a role in cellular host defense against Salmonella infection (4).
Synonyms:	LAMTOR2 Antibody, p14, ENDAP, ROBLD3, HSPC003, MAPBPIP, MAPKSP1AP, Ragulator2, Ragulator complex protein LAMTOR2, Endosomal adaptor protein p14, MAPBP-interacting protein
Host Species:	Rabbit
Clonality:	Polyclonal
Format:	IgG

Target Details

Gene Name:	LAMTOR2
Reactivity:	Human, Mouse, Rat
Immunogen Type:	Conjugated Peptide
Immunogen:	Anti-LAMTOR2 antibody was prepared from whole rabbit serum produced by repeated immunizations with an 18 amino acid peptide near the internal region of human LAMTOR2.

Purity/Specificity:	Anti-LAMTOR2 antibody was affinity purified from monospecific antiserum by immunoaffinity chromatography. LAMTOR2 antibody is human, mouse and rat reactive. At least two isoforms of LAMTOR2 are known to exist; this antibody will detect both isoforms. LAMTOR2 antibody is predicted to not cross-react with other LAMTOR family proteins.
Relevant Links:	• UniProtKB - Q9Y2Q5 • GeneID - 28956 • NCBI - NP_054736

Application Details

Tested Applications:	ELISA, IF, IHC, WB
Application Note:	Anti-LAMTOR2 Antibody has been tested for use in ELISA, Western Blotting, Immunohistochemistry and Immunofluorescence. Specific conditions for reactivity should be optimized by the end user. Expect a band at approximately 14 kDa in Western Blots of specific cell lysates and tissues.
Assay Dilutions:	All assays should be optimized by the user. Recommended dilutions (if any) may be listed below.
ELISA:	1:10,000 - 1:20,000
IF:	20 µg/mL
IHC:	5 µg/mL
WB:	1-2 µg/mL

Formulation

Physical State:	Liquid (sterile filtered)
Concentration:	1 mg/mL by UV absorbance at 280 nm
Buffer:	0.01 M Sodium Phosphate, 0.25 M Sodium Chloride, pH 7.2
Preservative:	0.02% (w/v) Sodium Azide
Stabilizer:	None


Shipping & Handling

Shipping Condition:	Dry Ice
----------------------------	---------

Storage Condition: Store vial at -20° C prior to opening. Aliquot contents and freeze at -20° C or below for extended storage. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use.

Expiration: Expiration date is one (1) year from date of receipt.

Images




Western Blot

Western blot of LAMTOR2.

Load: HepG2 cell lysate.

Primary Antibody: LAMTOR2 antibody at 1 µg/ml.




Immunohistochemistry

Immunohistochemistry of LAMTOR2.

Tissue: human spleen tissue.

Primary Antibody: LAMTOR2 antibody at 5 µg/mL.

**Immunofluorescence Microscopy**

Immunofluorescence of LAMTOR2.

Tissue: human spleen tissue.

Primary Antibody: LAMTOR2 antibody at 20 µg/mL.

Disclaimer

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 5199, Limerick, Pennsylvania, USA.