

Datasheet for 600-401-CV0**MMS21 Antibody****Overview**

Description:	Anti-MMS21 (RABBIT) Antibody - 600-401-CV0
Item No.:	600-401-CV0
Size:	100 µg
Applications:	ELISA, WB
Reactivity:	Human, Mouse, Rat
Host Species:	Rabbit

Product Details

Background:	MMS21, also known as NSE2, is a SUMO ligase that in combination with the SMC5/6 complex is required for the prevention of DNA damage induced apoptosis by facilitating DNA repair in human cells. MMS21-dependent sumoylation is integral and important to the cohesion mechanism and mitotic progression; this function appears to be independent of SMC6. MMS21 mediates SUMO attachment to various proteins such as SMC6L1 and TRAX, and possibly the cohesion components RAD21 and STAG2.
Synonyms:	MMS21 Antibody, NSE2, MMS21, ZMIZ7, C8orf36, E3 SUMO-protein ligase NSE2, MMS21 homolog, hMMS21
Host Species:	Rabbit
Clonality:	Polyclonal
Format:	IgG

Target Details

Gene Name:	NSMCE2
Reactivity:	Human, Mouse, Rat
Immunogen Type:	Conjugated Peptide
Immunogen:	Anti-MMS21 antibody was prepared from whole rabbit serum produced by repeated immunizations with a 15 amino acid synthetic peptide near the C-terminus of human MMS21.

Purity/Specificity:	Anti-MMS21 Antibody was affinity purified from monospecific antiserum by immunoaffinity chromatography. Cross reactivity with MMS21 from other sources has not been determined.
----------------------------	---

Relevant Links:	• UniProtKB - Q96MF7 • GeneID - 286053 • NCBI - NP_775956
------------------------	---

Application Details

Tested Applications:	ELISA, WB
-----------------------------	-----------

Application Note:	Anti-MMS21 Antibody has been tested for use in ELISA and Western Blotting. Specific conditions for reactivity should be optimized by the end user. Expect a band at approximately 28 kDa in Western Blots of specific cell lysates and tissues.
--------------------------	---

Assay Dilutions:	All assays should be optimized by the user. Recommended dilutions (if any) may be listed below.
-------------------------	---

ELISA:	1:5,000 - 1:10,000
---------------	--------------------

WB:	.5-1 µg/mL
------------	------------

Formulation

Physical State:	Liquid (sterile filtered)
------------------------	---------------------------

Concentration:	1 mg/mL by UV absorbance at 280 nm
-----------------------	------------------------------------

Buffer:	0.01 M Sodium Phosphate, 0.25 M Sodium Chloride, pH 7.2
----------------	---

Preservative:	0.02% (w/v) Sodium Azide
----------------------	--------------------------

Stabilizer:	None
--------------------	------


Shipping & Handling

Shipping Condition:	Dry Ice
----------------------------	---------

Storage Condition:	Store vial at -20° C prior to opening. Aliquot contents and freeze at -20° C or below for extended storage. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use.
---------------------------	---

Expiration:	Expiration date is one (1) year from date of receipt.
--------------------	---

Images



Western Blot

Western Blot of MMS21 antibody. Lane 1: 293 lysate with MMS21 antibody at 0.5 $\mu\text{g/mL}$. Lane 2: 293 lysate with MMS21 antibody at 1 $\mu\text{g/mL}$. Secondary antibody: Peroxidase rabbit secondary antibody at 1:10,000 for 45 min at RT. Block: 5% BLOTTO overnight at 4°C. Predicted/Observed size: 28 kDa, 30 kDa for MMS21. Other band(s): MMS21 splice variants and isoforms.

Disclaimer

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 5199, Limerick, Pennsylvania, USA.