

Datasheet for 600-401-D31**Dynamin Antibody****Overview**

Description:	Anti-Dynamin (RABBIT) Antibody - 600-401-D31
Item No.:	600-401-D31
Size:	100 µL
Applications:	WB
Reactivity:	Rat
Host Species:	Rabbit

Product Details

Background:	Dynamin Antibody detects Dynamin which is a member of a group of nerve terminal proteins called dephosphins that regulate synaptic vesicle endocytosis. There are 3 known isoforms of Dynamin, each having several splice variants as well. Dynamin I is expressed only in neurons whereas Dynamin II is ubiquitously expressed and Dynamin III is found primarily in the testes. Dynamin 1 is phosphorylated by PKC and dephosphorylated by calcineurin. Anti-Dynamin Antibody is ideal for investigators involved in Cell Signaling, Cell Motility, Neuroscience and Cancer Research.
Synonyms:	Dynamin-1, B-dynamin, D100
Host Species:	Rabbit
Clonality:	Polyclonal
Format:	IgG

Target Details

Gene Name:	DNM1
Reactivity:	Rat
Immunogen Type:	Conjugated Peptide
Immunogen:	Anti-Dynamin Antibody was produced in rabbit by repeated immunizations with synthetic peptide corresponding to amino acid residues from the N-terminal region conjugated to KLH.

Purity/Specificity: Anti-Dynamin antibody is directed against human Dynamin. Anti-Dynamin antibody was affinity purified from monospecific antiserum by immunoaffinity purification. Anti-Dynamin antibody recognizes Dynamin from rat sources. Cross reactivity with Dynamin from other species has not been determined. However, reactivity is also expected from the following species based on 100% sequence homology: bovine, chicken, human, mouse and non-human primate.

Relevant Links:

- [UniProtKB - Q05193](#)
- [GeneID - 1759](#)
- [NCBI - AAA88025.1](#)

Application Details

Tested Applications: WB

Application Note: Anti-Dynamin Antibody is tested for use in ELISA and Western Blotting. Specific conditions for reactivity should be optimized by the end user. Expect a band of approximately 95 kDa in size corresponding to the dynamin proteins in the appropriate cell lysate or extract.

Assay Dilutions: All assays should be optimized by the user. Recommended dilutions (if any) may be listed below.

ELISA: 1:10,000

WB: 1:1000

Formulation

Physical State: Liquid

Buffer: 0.01 M HEPES, 0.15 M Sodium Chloride, pH 7.5

Stabilizer: 0.1 mg/ml Bovine Serum Albumin (BSA) - IgG and Protease free, 50% (v/v) Glycerol


Shipping & Handling

Shipping Condition: Dry Ice

Storage Condition: Store vial at -20° C prior to opening. This product is stable at 4° C as an undiluted liquid. For extended storage, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Dilute only prior to immediate use.

Expiration: Expiration date is one (1) year from date of receipt.

Images



Western Blot

Western Blot of Anti-Dynamin (Rabbit) Antibody. Lane 1: rat cortex lysate. Lane 2: none. Load: 10 µg per lane. Primary antibody: Dynamin antibody at 1:400 for overnight at 4°C. Secondary antibody: IRDye800™ rabbit secondary antibody at 1:10,000 for 45 min at RT. Block: 5% BLOTTO overnight at 4°C. Predicted/Observed size: ~95kDa/~95kDa for dynamin protein. Other band(s): none.

Disclaimer

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 5199, Limerick, Pennsylvania, USA.