

Datasheet for 600-401-DE6**N-RAS Antibody****Overview**

Description:	Anti-N-RAS (RABBIT) Antibody - 600-401-DE6
Item No.:	600-401-DE6
Size:	100 µg
Applications:	ELISA, WB
Reactivity:	Human, Mouse, Rat
Host Species:	Rabbit

Product Details

Background:	Activating mutations and overexpression of classical Ras subfamily members (N-RAS, H-RAS, and K-RAS) have been widely investigated as key events in the development of human cancers. The N-RAS protein shuttles between the Golgi apparatus and the plasma membrane. This shuttling is regulated through palmitoylation and depalmitoylation by the ZDHHC9-GOLGA7 complex. N-RAS, which has intrinsic GTPase activity, is activated by a guanine nucleotide-exchange factor and inactivated by a GTPase activating protein. Mutations in this gene have been associated with somatic rectal cancer, follicular thyroid cancer, autoimmune lymphoproliferative syndrome, Noonan syndrome, and juvenile myelomonocytic leukemia.
Synonyms:	N-RAS Antibody, NS6, ALPS4, N-ras, NRAS1, HRAS1, GTPase NRas, Transforming protein N-Ras
Host Species:	Rabbit
Clonality:	Polyclonal
Format:	IgG

Target Details

Gene Name:	NRAS
Reactivity:	Human, Mouse, Rat
Immunogen Type:	Conjugated Peptide
Immunogen:	Anti-N-RAS antibody was prepared from whole rabbit serum produced by repeated immunizations with a 15 amino acid peptide near the C-terminus of human N-RAS.

Purity/Specificity:	Anti-N-RAS Antibody was affinity purified from monospecific antiserum by immunoaffinity chromatography. N-RAS antibody is predicted to not cross-react with H-RAS
----------------------------	---

Relevant Links:	• UniProtKB - P01111 • GeneID - 4893 • NCBI - NP_002515
------------------------	---

Application Details

Tested Applications:	ELISA, WB
-----------------------------	-----------

Application Note:	Anti-N-RAS Antibody has been tested for use in ELISA and Western Blotting. Expect a band at approximately 21 kDa in Western Blots using positive control A431 whole cell lysates or other specific cell lysates and tissues. Specific conditions for reactivity should be optimized by the end user.
--------------------------	--

Assay Dilutions:	All assays should be optimized by the user. Recommended dilutions (if any) may be listed below.
-------------------------	---

ELISA:	1:10,000 - 1:20,000
---------------	---------------------

WB:	1 and 2 µg/mL
------------	---------------

Formulation

Physical State:	Liquid (sterile filtered)
------------------------	---------------------------

Concentration:	1 mg/mL by UV absorbance at 280 nm
-----------------------	------------------------------------

Buffer:	0.01 M Sodium Phosphate, 0.25 M Sodium Chloride, pH 7.2
----------------	---

Preservative:	0.02% (w/v) Sodium Azide
----------------------	--------------------------

Stabilizer:	None
--------------------	------


Shipping & Handling

Shipping Condition:	Wet Ice
----------------------------	---------

Storage Condition:	Antibody can be stored at 4°C up to one year. Antibodies should not be exposed to prolonged high temperatures.
---------------------------	--

Expiration:	Expiration date is one (1) year from date of receipt.
--------------------	---

Images



Western Blot

Western blot of NRAS.

Load: A431 cell lysate.

Primary Antibody: NRAS antibody at (A) 1 and (B) 2 $\mu\text{g/mL}$.

Disclaimer

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 5199, Limerick, Pennsylvania, USA.