

Datasheet for 600-401-EH3**SFTS Virus HB29 Antibody****Overview**

Description:	Anti-SFTS Virus HB29 (RABBIT) Antibody - 600-401-EH3
Item No.:	600-401-EH3
Size:	100 µg
Applications:	ELISA
Reactivity:	Virus
Host Species:	Rabbit

Product Details

Background:	Severe fever with thrombocytopenia syndrome (SFTS) is an emerging infectious disease caused by SFTS virus, a newly discovered bunyavirus that can cause high rates of fatalities. This disease is thought to be transmitted through blood contact. The SFTS virus HB29 membrane glycoprotein polypeptide mRNA encodes two glycoproteins termed Phlebovirus glycoprotein G1 and G2 respectively. This antibody will detect the non-cleaved glycoprotein
Synonyms:	SFTS Virus HB29 Antibody
Host Species:	Rabbit
Clonality:	Polyclonal
Format:	IgG

Target Details

Gene Name:	B303_sMgp1
Reactivity:	Virus
Immunogen Type:	Conjugated Peptide
Immunogen:	Anti-SFTS Virus HB29 Membrane Glycoprotein antibody was prepared from whole rabbit serum produced by repeated immunizations with a 15 amino acid synthetic peptide near the C-terminus of human SFTS Virus HB29 Membrane Glycoprotein.

Purity/Specificity:	Anti-SFTS Virus HB29 Antibody was affinity purified from monospecific antiserum by immunoaffinity chromatography. SFTS Virus HB29 Membrane Glycoprotein antibody is specific to the SFTS Virus HB29.
Relevant Links:	• UniProtKB - F1BA47 • GeneID - 13231111 • NCBI - ADZ04471

Application Details

Tested Applications:	ELISA
Application Note:	Anti-SFTS Virus HB29 Antibody has been tested for use in ELISA. Specific conditions for reactivity should be optimized by the end user. Expect a band at approximately 117 kDa in Western Blots of specific cell lysates and tissues.
Assay Dilutions:	All assays should be optimized by the user. Recommended dilutions (if any) may be listed below.
ELISA:	User Optimized

Formulation

Physical State:	Liquid (sterile filtered)
Concentration:	1.0 mg/mL by UV absorbance at 280 nm
Buffer:	0.01 M Sodium Phosphate, 0.25 M Sodium Chloride, pH 7.2
Preservative:	0.02% (w/v) Sodium Azide
Stabilizer:	None

Shipping & Handling

Shipping Condition:	Wet Ice
Storage Condition:	Antibody can be stored at 4°C up to one year. Antibodies should not be exposed to prolonged high temperatures.
Expiration:	Expiration date is one (1) year from date of receipt.

Disclaimer

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 5199, Limerick, Pennsylvania, USA.