

Datasheet for 600-401-EK8**Sin Nombre Virus Glycoprotein 2 Antibody****Overview**

Description:	Anti-Sin Nombre Virus Glycoprotein 2 (RABBIT) Antibody - 600-401-EK8
Item No.:	600-401-EK8
Size:	100 µg
Applications:	ELISA
Reactivity:	Virus
Host Species:	Rabbit

Product Details

Background:	Sin Nombre virus (SNV) is a rodent-borne hantavirus of the family Bunyaviridae, an enveloped, negative-sense RNA viruses with a tripartite genome that can cause hantavirus pulmonary syndrome (HPS) (1). Hantavirus glycoprotein precursor (GPC) is posttranslationally cleaved into two glycoproteins G1 (Gn) and G2 (Gc). While the G1 glycoprotein is thought to be degraded by the host autophagy machinery, and this autophagic clearance is required for efficient virus replication (2), no such degradation was observed for the G2 glycoprotein.
Synonyms:	Sin Nombre Virus Glycoprotein 2 Antibody
Host Species:	Rabbit
Clonality:	Polyclonal
Format:	IgG

Target Details

Gene Name:	SNVsMgp1
Reactivity:	Virus
Immunogen Type:	Conjugated Peptide
Immunogen:	Anti-Sin Nombre virus glycoprotein 2 antibody was prepared from whole rabbit serum produced by repeated immunizations with a 16 amino acid peptide near the internal region of the Sin Nombre virus glycoprotein.

Purity/Specificity: Anti-Sin Nombre virus glycoprotein 2 antibody was affinity purified from monospecific antiserum by immunoaffinity chromatography. Cross reactivity with Sin Nombre Virus Glycoprotein 2 from other sources has not been determined.

Relevant Links:

- [UniProtKB - Q89905](#)
- [GeneID - 2654026](#)
- [NCBI - NP_941974](#)

Application Details

Tested Applications: ELISA

Application Note: Anti-Sin Nombre Virus Glycoprotein 2 Antibody has been tested for use in ELISA. Specific conditions for reactivity should be optimized by the end user. Expect a band at approximately 126 kDa in Western Blots of specific cell lysates and tissues.

Assay Dilutions: All assays should be optimized by the user. Recommended dilutions (if any) may be listed below.

ELISA: User Optimized

Formulation

Physical State: Liquid (sterile filtered)

Concentration: 1 mg/mL by UV absorbance at 280 nm

Buffer: 0.01 M Sodium Phosphate, 0.25 M Sodium Chloride, pH 7.2

Preservative: 0.02% (w/v) Sodium Azide

Stabilizer: None

Shipping & Handling

Shipping Condition: Dry Ice

Storage Condition: Store vial at -20° C prior to opening. Aliquot contents and freeze at -20° C or below for extended storage. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use.

Expiration: Expiration date is one (1) year from date of receipt.

Disclaimer

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 5199, Limerick, Pennsylvania, USA.