

Datasheet for 600-401-FV9**VARP Antibody****Overview**

Description:	Anti-VARP (RABBIT) Antibody - 600-401-FV9
Item No.:	600-401-FV9
Size:	100 µg
Applications:	ELISA, IHC
Reactivity:	Human
Host Species:	Rabbit

Product Details

Background:	The VPS9 ankyrin repeat protein (VARP) binds to the Rab21, a guanine nucleotide exchange factor that plays an essential role in endocytic trafficking. VARP localizes to early endosomes and is thought to regulate endosome dynamics. VARP also interacts with TI-VAMP/VAMP7, a vesicular SNARE that mediates an exocytic pathway that is crucial to neurite growth. Depletion of VARP by RNA interference impairs neurite growth, suggesting that VARP is a positive regulator of neurite growth.
Synonyms:	VARP Antibody, VARP, PP12899, Ankyrin repeat domain-containing protein 27, VPS9 domain-containing protein
Host Species:	Rabbit
Clonality:	Polyclonal
Format:	IgG

Target Details

Gene Name:	ANKRD27
Reactivity:	Human
Immunogen Type:	Conjugated Peptide
Immunogen:	Anti-VARP antibody was prepared from whole rabbit serum produced by repeated immunizations with a 15 amino acid synthetic peptide near the internal region of human VARP.

Purity/Specificity: Anti-VARP Antibody was affinity purified from monospecific antiserum by immunoaffinity chromatography. Cross reactivity with VARP from other sources has not been determined.

Relevant Links:

- [UniProtKB - Q96NW4](#)
- [GeneID - 84079](#)
- [NCBI - NP_115515.2](#)

Application Details

Tested Applications: ELISA, IHC

Application Note: Anti-VARP Antibody has been tested for use in ELISA and Immunohistochemistry. Specific conditions for reactivity should be optimized by the end user. Expect a band at approximately 117 kDa in Western Blots of specific cell lysates and tissues.

Assay Dilutions: All assays should be optimized by the user. Recommended dilutions (if any) may be listed below.

ELISA: User Optimized

IHC: 5 µg/mL

Formulation

Physical State: Liquid (sterile filtered)

Concentration: 1 mg/mL by UV absorbance at 280 nm

Buffer: 0.01 M Sodium Phosphate, 0.25 M Sodium Chloride, pH 7.2

Preservative: 0.02% (w/v) Sodium Azide

Stabilizer: None


Shipping & Handling

Shipping Condition: Dry Ice

Storage Condition: Store vial at -20° C prior to opening. Aliquot contents and freeze at -20° C or below for extended storage. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use.

Expiration: Expiration date is one (1) year from date of receipt.

Images



Immunohistochemistry

Immunohistochemistry of VARP antibody. Tissue: Human small intestine. Fixation: formalin fixed paraffin embedded. Antigen retrieval: not required. Primary antibody: VARP antibody at 5 µg/mL for 1 h at RT. Secondary antibody: Peroxidase rabbit secondary antibody at 1:10,000 for 45 min at RT. Localization: VARP is a multi-pass membrane protein. Staining: VARP as precipitated pink signal with blue nuclear counterstain.

Disclaimer

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 5199, Limerick, Pennsylvania, USA.