

Datasheet for 600-401-H23**Swine H1N1 Hemagglutinin Antibody****Overview**

Description:	Anti-Swine H1N1 Hemagglutinin (RABBIT) Antibody - 600-401-H23
Item No.:	600-401-H23
Size:	100 µg
Applications:	ELISA
Reactivity:	Virus
Host Species:	Rabbit

Product Details**Background:**

Swine H1N1 Hemagglutinin Antibody detects novel swine influenza Hemagglutinin protein. Influenza A virus is a major public health threat, killing more than 30,000 people per year in the USA. In early 2009, a novel swine-origin influenza A (H1N1) virus was identified in specimens obtained from patients in Mexico and the United States. The virus spread quickly around the world and on June 11, 2009, the World Health Organization declared it a pandemic. Influenza A virus has one of sixteen possible Hemagglutinin (HA) surface proteins and one of nine possible Neuraminidase (NA) surface proteins. The Hemagglutinin protein facilitates viral attachment while Neuraminidase is involved in viral release. These proteins also elicit immune responses that prevent infection or independently reduce viral replication. The genetic make-up of this swine flu virus is unlike any other: it is an H1N1 strain that combines a triple assortment first identified in 1998 including human, swine, and avian influenza with two new pig H3N2 virus genes from Eurasia, themselves of recent human origin. The distinct antigenic properties of the new swine virus compared with seasonal influenza A (H1N1) virus suggest that human immunity against new swine influenza virus is limited, although the age distribution of reported cases suggests some degree of protection in older age groups. Anti-Swine H1N1 Hemagglutinin antibodies are ideal for investigators involved in infectious disease research.

Synonyms:	Swine-Origin Influenza A (H1N1) Hemagglutinin Antibody, Swine-Origin Influenza A (H1N1) Hemagglutinin, S-OIV A Hemagglutinin, Swine flu H1, HA
Host Species:	Rabbit
Clonality:	Polyclonal
Format:	IgG

Target Details

Gene Name:	HA
Reactivity:	Virus
Immunogen Type:	Conjugated Peptide
Immunogen:	Swine H1N1 Hemmagglutinin Antibody was produced from whole rabbit serum prepared by repeated immunizations with a synthetic peptide from the novel swine influenza Hemagglutinin protein.
Purity/Specificity:	Anti-Swine H1N1 Hemmagglutinin Antibody was affinity purified from monospecific antiserum by immunoaffinity purification. A BLAST analysis was used to suggest cross-reactivity with Swine H1N1 Hemagglutinin with virus based on 100% homology with the immunizing sequence. Cross-reactivity with Swine H1N1 Hemagglutinin from other sources has not been determined. This antibody is specific for the novel Swine influenza Hemagglutinin and will not recognize the corresponding Hemagglutinin sequence from seasonal H1N1 influenza (A/Brisbane/49/2007 (H1N1)).
Relevant Links:	• NCBI - ACQ76314 • UniProtKB - C4AL34

Application Details

Tested Applications:	ELISA
Application Note:	Anti-Swine H1N1 Hemmagglutinin A Antibody is tested for use in ELISA and WB. Expect a band approximately ~63.2 kDa on specific lysates. Specific conditions for reactivity should be optimized by the end user.
Assay Dilutions:	All assays should be optimized by the user. Recommended dilutions (if any) may be listed below.
ELISA:	1 µg/ml


Formulation

Physical State:	Liquid (sterile filtered)
Concentration:	1 mg/mL by UV absorbance at 280 nm
Buffer:	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2
Preservative:	0.02% (w/v) Sodium Azide
Stabilizer:	None

Shipping & Handling

Shipping Condition:	Dry Ice
Storage Condition:	Store vial at -20° C prior to opening. Aliquot contents and freeze at -20° C or below for extended storage. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use.
Expiration:	Expiration date is one (1) year from date of receipt.

Images



Western Blot

Western blot of Anti-Swine H1N1 Hemagglutinin.

Load: recombinant seasonal hemagglutinin (lanes A & B) and swine-origin hemagglutinin (lanes C & D).

Primary Antibody: Anti-swine H1N1 Hemagglutinin antibody (p/n 600-401-H23) at 2 µg/mL (lanes B & D) and a anti-seasonal Hemagglutinin antibody at 2 µg/mL (lanes A & C).

Disclaimer

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 5199, Limerick, Pennsylvania, USA.