

Datasheet for 600-401-J52**Tropomyosin-1 alpha phospho S283 Antibody****Overview**

Description:	Anti-Tropomyosin-1 alpha phosphoS283 (RABBIT) Antibody - 600-401-J52
Item No.:	600-401-J52
Size:	50 µg
Reactivity:	Human
Host Species:	Rabbit

Product Details

Background:	Anti-Tropomyosin-1 alpha pS283 antibody attaches to actin filaments in muscle and non-muscle cells, a large factor in the calcium dependent regulation of muscle contraction. Anti-Tropomyosin-1 alpha pS283 antibody is ideal for researchers interested in Signal Transduction and Cardiovascular research.
Synonyms:	rabbit anti-Tropomyosin-1 alpha pS283 antibody, rabbit anti-Tropomyosin1 alpha pS283 antibody, Tropomyosin α -1 chain, Alpha-tropomyosin, α -tropomyosin, Tropomyosin-1, TPM1, TPM-1, C15orf13, TMSA
Host Species:	Rabbit
Clonality:	Polyclonal
Format:	IgG

Target Details

Gene Name:	TPM1
Reactivity:	Human
PTM Specificity:	Phosphorylation
Immunogen Type:	Conjugated Peptide
Immunogen:	Tropomyosin-1 alpha pS283 affinity purified antibody was prepared from whole rabbit serum produced by repeated immunizations with a synthetic phospho-peptide surrounding S283 human Tropomyosin-1 alpha pS283.

Purity/Specificity: Anti-Tropomyosin-1 alpha pS283 was affinity purified from monospecific antiserum by immunoaffinity chromatography and is directed against the phosphorylated form of human S283 residue. A BLAST analysis was used to suggest cross-reactivity with mouse and human based on 100% sequence homology. Cross-reactivity with Tropomyosin-1 alpha pS283 from other sources has not been determined.

Relevant Links:

- [UniProtKB - P09493](#)
- [NCBI - NP_000357.3](#)
- [GeneID - 7168](#)

Application Details

Application Note: Anti-Tropomyosin-1 alpha pS283 antibody is useful for ELISA, immunohistochemistry, and Western Blot. Specific conditions for reactivity should be optimized by the end user. Expect a band approximately ~32.7kDa corresponding to the appropriate cell lysate or extract.

Assay Dilutions: All assays should be optimized by the user. Recommended dilutions (if any) may be listed below.

ELISA: 1: 10,000

WB: 1ug/ml

Formulation

Physical State: Lyophilized

Concentration: 0.80mg/mL by UV absorbance at 280 nm

Buffer: 0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2

Stabilizer: 10 mg/mL Bovine Serum Albumin (BSA) - Immunoglobulin and Protease free

Reconstitution Volume: 50μL

Reconstitution Buffer: Restore with deionized water (or equivalent)

Shipping & Handling

Shipping Condition: Ambient

Storage Condition: Store vial at 4° C prior to restoration. For extended storage aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use.

Expiration: Expiration date is one (1) year from date of receipt.

Disclaimer

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 5199, Limerick, Pennsylvania, USA.