

Datasheet for 600-401-J80**Huntington phospho S1181 Antibody****Overview**

Description:	Anti-Huntington pS1181 (RABBIT) Antibody - 600-401-J80
Item No.:	600-401-J80
Size:	50 µg
Reactivity:	Human, Mouse, Rat
Host Species:	Rabbit

Product Details

Background:	Huntington's disease (HD) is an autosomal dominant neurological disorder caused by a polyglutamine (polyQ) repeat expansion in the huntingtin (Htt) protein. The disease is characterized by neurodegeneration and formation of neuronal intracellular inclusions primarily in the striatum and cortex, leading to personality changes, motor impairment, and dementia. The Huntington's disease protein is ~350kD in size and is localized in the brain. Anti-Huntington antibody is ideal for investigators interested in Neuroscience research.
Synonyms:	rabbit anti-Huntington pS1181 antibody, rabbit anti-Huntingtin pS1181 antibody, HD protein, HTT, Huntington Disease Protein, Huntington disease, IT15, LOMARS
Host Species:	Rabbit
Clonality:	Polyclonal
Format:	IgG

Target Details

Gene Name:	HTT
Reactivity:	Human, Mouse, Rat
PTM Specificity:	Phosphorylation
Immunogen Type:	Conjugated Peptide
Immunogen:	Anti-Huntingtin pS1181 affinity purified antibody was prepared from whole rabbit serum produced by repeated immunizations with a synthetic phospho-peptide corresponding to the internal region of human Huntington disease protein.

Purity/Specificity: Anti-Huntingtin pS1181 was affinity purified from monospecific antiserum by immunoaffinity chromatography and is directed against the phosphorylated form of the S1181 residue. A BLAST analysis was used to suggest cross-reactivity with human, mouse, and rat based on 100% sequence homology. Cross-reactivity with Huntingtin pS1181 from other sources has not been determined.

Relevant Links:

- [NCBI - NP_002102.4](#)
- [GeneID - 3064](#)
- [UniProtKB - P42858](#)

Application Details

Application Note: Anti-HTT antibody is useful for ELISA, Immunohistochemistry, and Western Blot. Specific conditions for reactivity should be optimized by the end user. Expect a band approximately ~350kDa corresponding to the appropriate cell lysate or extract.

Assay Dilutions: All assays should be optimized by the user. Recommended dilutions (if any) may be listed below.

ELISA: 1: 10,000

IHC: User Optimized

WB: 1ug/ml

Formulation

Physical State: Liquid (sterile filtered)

Concentration: 1.0 mg/ml by UV absorbance at 280 nm

Buffer: 0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2

Preservative: None

Stabilizer: 50% (v/v) Glycerol with 1 mg/ml Bovine Serum Albumin (BSA)

Shipping & Handling

Shipping Condition: Dry Ice

Storage Condition: Store vial at -20° C prior to opening. Aliquot contents and freeze at -20° C or below for extended storage. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use.

Expiration: Expiration date is one (1) year from date of receipt.

Disclaimer

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 5199, Limerick, Pennsylvania, USA.