

Datasheet for 600-401-MG0S**Smoothened Antibody****Overview**

Description:	Anti-Smoothened (RABBIT) Antibody - 600-401-MG0S
Item No.:	600-401-MG0S
Size:	25 µL
Applications:	ELISA, IHC, WB
Reactivity:	Human
Host Species:	Rabbit

Product Details

Background:	SMO (Smoothened, Frizzled Class Receptor) receptors are non-classical G protein-coupled receptors that probably associates with the patched protein (PTCH) to transduce the hedgehogs proteins signal. The binding of sonic hedgehog (SHH) to its receptor patched is thought to prevent normal inhibition by patched of smoothened (SMO). It is required for the accumulation of KIF7, GLI2 and GLI3 in the cilia. And it interacts with DLG5 at the ciliary base to induce the accumulation of KIF7 and GLI2 at the ciliary tip for GLI2 activation. Diseases associated with SMO include Curry-Jones Syndrome and Basal Cell Carcinoma. Anti-Smoothened Antibody is useful for researchers interested in Cancer research.
Synonyms:	Rabbit Anti-SMO Antibody, Smoothened Frizzled Class Receptor, Smoothened Seven Transmembrane Spanning Receptor, Smoothened Frizzled Family Receptor, Frizzled Family Member 11, Protein GX, SMOH, Seven Transmembrane Helix Receptor, Smoothened (Drosophila) Homolog, Smoothened Homolog (Drosophila), Smoothened Homolog, FZD11, CRJS, Gx
Host Species:	Rabbit
Clonality:	Polyclonal
Format:	IgG

Target Details

Gene Name:	SMO
Reactivity:	Human

Immunogen Type:	Conjugated Peptide
Immunogen:	Anti-Smoothened antibody was prepared from whole rabbit serum produced by repeated immunizations with a synthetic peptide corresponding to an internal portion of human Smoothened conjugated to Keyhole Limpet Hemocyanin (KLH).
Purity/Specificity:	This affinity purified antibody is directed against human SMO. This product was affinity purified from monospecific antiserum by immunoaffinity purification. A BLAST analysis was used to suggest cross-reactivity with the antigen based on 100% homology with the immunizing sequence to mouse and rat.
Relevant Links:	• UniProtKB - Q99835 • NCBI - NP_005622.1 • GeneID - 6608

Application Details

Tested Applications:	ELISA, IHC, WB
Application Note:	Anti-Smoothened Antibody has been tested in ELISA, WB, and IHC. Expect a band at ~86.4kDa in western blot using appropriate lysates and tissues. Positive control used: Overexpressing lysate and HEK293T in Western Blot; and Human Lung Tissue in Immunohistochemistry.
Assay Dilutions:	All assays should be optimized by the user. Recommended dilutions (if any) may be listed below.
ELISA:	1:10,000 - 1:50,000
IHC:	1:100
WB:	1 µg/ml

Formulation

Physical State:	Liquid (sterile filtered)
Concentration:	1.0 mg/ml by UV absorbance at 280 nm
Buffer:	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2
Preservative:	0.01% (w/v) Sodium Azide
Stabilizer:	None


Shipping & Handling

Shipping Condition:	Dry Ice
----------------------------	---------



Storage Condition: Store vial at -20° C or below prior to opening. This vial contains a relatively low volume of reagent (25 µL). To minimize loss of volume dilute 1:10 by adding 225 µL of the buffer stated above directly to the vial. Recap, mix thoroughly and briefly centrifuge to collect the volume at the bottom of the vial. Use this intermediate dilution when calculating final dilutions as recommended below. Store the vial at -20°C or below after dilution. Avoid cycles of freezing and thawing.

Expiration: Expiration date is one (1) year from date of receipt.

Images




A. Human lung stained with H & E at 200x
 B. Human lung stained with Anti-SMO at 200x, 1:200 dilution
 C. Human lung stained with Anti-SMO at 400x, 1:200 dilution



Immunohistochemistry

Immunohistochemistry of Rabbit Anti-Smoothed Antibody. Tissue: Human Lung Tissue. Fixative: None. Antigen Retrieval: HIER using Citrate buffer for 20 mins. Staining: (A) H&E (200x), (B) Anti-SMO at 1:200 (200x), (C) Anti-SMO at 1:200 (400x). Secondary Antibody: Anti-Rabbit Poly IgG HRP for 8mins at RT. Counterstain: Hematoxylin. Results: SMO shows positive nuclear staining in alveolar cells and bronchial epithelium cells.



Western Blot

Western Blot of Rabbit Anti-Smoothed Antibody. Lane 1: Opal Prestained Molecular Weight Marker (p/n MB-210-0500). Lane 2: Smoothed Overexpressed HEK293T Lysate DDK-Tagged [+]. Lane 3: HEK293T Lysate (p/n W09-001-GX5). Load: 10µg/lane. Primary Antibody: Anti-SMO at 1:1000 overnight at 2-8°C. Secondary Antibody: Goat anti-Rabbit IgG HRP (p/n 611-103-122) 1:70,000 for 1hr at RT. Blocking Buffer: Universal BlockOut (p/n MB-073). Predicted MW: ~100kDa overexpressed lysate. Observed MW: ~48, 110kDa. Exposure: 9sec.

Disclaimer

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 5199, Limerick, Pennsylvania, USA.