

Datasheet for 600-401-MG1**SLC40A1 Antibody****Overview**

Description:	Anti-SLC40A1 (RABBIT) Antibody - 600-401-MG1
Item No.:	600-401-MG1
Size:	100 µg
Applications:	ELISA, IF, WB
Reactivity:	Human
Host Species:	Rabbit

Product Details

Background:	SLC40A1 (Solute Carrier Family 40 Member 1) may be involved in iron export from duodenal epithelial cell and in transfer of iron between maternal and fetal circulation. Mediates iron efflux in the presence of a ferroxidase (hephaestin and/or ceruloplasmin). Defects in this gene are a cause of hemochromatosis type 4 (HFE4). Diseases associated with SLC40A1 include Hemochromatosis Type 4 and Type 1. Anti-SLC40A1 Antibody is useful for researchers interested in mineral absorption, insulin receptor recycling, and transmembrane transporter activity.
Synonyms:	Rabbit Anti-SLC40A1 Antibody, Solute Carrier Family 40 Member 1, Solute Carrier Family 11 (Proton-Coupled Divalent Metal Ion Transporters), Member 3, Solute Carrier Family 40 (Iron-Regulated Transporter), Member 1, SLC11A3, IREG1, FPN1, Iron-Regulated Transporter 1, Iron Regulated Gene 1, Ferroportin-1, Ferroportin 1, MSTP079, MST079, HFE4, MTP1
Host Species:	Rabbit
Clonality:	Polyclonal
Format:	IgG

Target Details

Gene Name:	SLC40A1
Reactivity:	Human
Immunogen Type:	Conjugated Peptide

Immunogen:	Anti-SLC40A1 antibody was prepared from whole rabbit serum produced by repeated immunizations with a synthetic peptide corresponding to an internal portion of human SLC40A1 conjugated to Keyhole Limpet Hemocyanin (KLH).
Purity/Specificity:	This affinity purified antibody is directed against human SLC40A1. This product was affinity purified from monospecific antiserum by immunoaffinity purification. A BLAST analysis was used to suggest cross-reactivity with the antigen based on 100% homology with the immunizing sequence to mouse and rat.
Relevant Links:	• UniProtKB - Q9NP59 • NCBI - NP_055400.1 • GeneID - 30061

Application Details

Tested Applications:	ELISA, IF, WB
Application Note:	Anti-SLC40A1 Antibody has been tested in ELISA, WB, and IF. Expect a band at ~62.5kDa in western blot using appropriate lysates. Positive control used: Caco-2 lysate in Western Blot and Hep-G2 in Immunofluorescence.
Assay Dilutions:	All assays should be optimized by the user. Recommended dilutions (if any) may be listed below.
ELISA:	1:10,000 - 1:50,000
IF:	15 µg/ml
WB:	10 µg/ml

Formulation

Physical State:	Liquid (sterile filtered)
Concentration:	1.04 mg/ml by UV absorbance at 280 nm
Buffer:	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2
Preservative:	0.01% (w/v) Sodium Azide
Stabilizer:	None


Shipping & Handling

Shipping Condition:	Dry Ice
----------------------------	---------

Storage Condition: Store vial at -20° C prior to opening. Aliquot contents and freeze at -20° C or below for extended storage. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use.


Expiration: Expiration date is one (1) year from date of receipt.

Images




Immunofluorescence Microscopy

Immunofluorescence of Rabbit Anti-SLC40A1 Antibody. Cells: Hep-G2 Cells. Fixative: 4% PFA. Permeabilization: 0.3% Triton X-100. Primary Antibody: Anti-SLC40A1 at 15μg/mL overnight at 2-8°C. Counterstain: Donkey Anti-Rabbit IgG DyLight™488 (p/n 611-741-127) at 5μg/mL for 1hr at RT. Nuclear Counterstain: DAPI. Staining: (A) DAPI. (B) Anti-SLC40A1 + DyLight™488. (C) Merge A + B. (D) secondary only. Expected Localization: cytoplasm/plasma.



Western Blot

Western Blot of Rabbit Anti-SLC40A1 Antibody. Lane 1: Opal Pre-stained molecular weight (p/n MB-210-0500). Lane 2: Caco-2 Lysate (p/n W09-001-GX4). Lane 3: K-562 Lysate (p/n W09-001-GJ7). Lane 4: Hep-G2 Lysate (p/n W09-001-GJ5). Lane 5: MCF-7 Lysate (p/n W09-000-360) [-]. Primary Antibody: Anti-SLC40A1 at 5μg/mL overnight at 2-8°C. Secondary Antibody: Goat anti-Rabbit IgG HRP (p/n 611-103-122) at 1:70,000 for 1 hr at RT. Blocking Buffer: Universal BlockOut Buffer (p/n MB-073). Predicted MW: ~63kDa. Observed MW: ~75kDa. Exposure time: 40 seconds.



Western Blot

Western Blot of Rabbit Anti-SLC40A1 Antibody. Lane 1: Opal Pre-stained molecular weight (p/n MB-210-0500). Lane 2: SLC40A1 overexpressing HEK293T Lysate (10µg). Lane 3: HEK293T empty vector Lysate (10µg). Primary Antibody: Anti-SLC40A1 at 1:1000 overnight at 2-8°C. Secondary Antibody: Goat anti-Rabbit IgG HRP (p/n 611-103-122) at 1:70,000 for 1 hr at RT. Blocking Buffer: Universal BlockOut Buffer (p/n MB-073). Predicted MW: ~70kDa for overexpressed. Observed MW: ~75kDa. Exposure time: 38 seconds.

Disclaimer

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 5199, Limerick, Pennsylvania, USA.