

Datasheet for 600-401-P10**Histone H3 K79me1 Antibody****Overview**

Description:	Anti-Histone H3 [Monomethyl Lys79] (RABBIT) Antibody - 600-401-P10
Item No.:	600-401-P10
Size:	50 µg
Applications:	ChIP, Dot Blot, IF, Multiplex, WB
Reactivity:	Human, Mouse, Monkey
Host Species:	Rabbit

Product Details

Background:	Chromatin is the arrangement of DNA and proteins in which chromosomes are formed. Correspondingly, chromatin is formed from nucleosomes, which are comprised of a set of four histone proteins (H2A, H2B, H3, H4) wrapped with DNA. Chromatin is a very dynamic structure in which numerous post-translational modifications work together to activate or repress the availability of DNA to be copied, transcribed, or repaired. These marks decide which DNA will be open and commonly active (euchromatin) or tightly wound to prevent access and activation (heterochromatin). Common histone modifications include methylation of lysine and arginine, acetylation of lysine, phosphorylation of threonine and serine, and sumoylation, biotinylation, and ubiquitylation of lysine. In particular Lys79 methylations are involved in transcriptional activation, and has been found to be inversely correlated with H2B ubiquitination. Anti-Histone H3 K79-Me1 antibody is ideal for researchers interested in Chromatin Modifiers, Chromatin Research, Histones and Modified Histones, and Epigenetics research.
Synonyms:	rabbit anti-Histone H3 monomethyl Lys79 antibody, H3K79me1, Anti-H3 K79-Me1, Histone H3 antibodies
Host Species:	Rabbit
Clonality:	Polyclonal
Format:	IgG

Target Details

Gene Name:	HIST1H3A
Reactivity:	Human, Mouse, Monkey

PTM Specificity:	Methylation
Immunogen Type:	Conjugated Peptide
Immunogen:	Histone H3 K79-Me1 affinity purified antibody was prepared from whole rabbit serum produced by repeated immunizations with a synthetic peptide surrounding the Lys79 site of human Histone H3.
Purity/Specificity:	Anti-Histone H3 [Monomethyl Lys79] was affinity purified from monospecific antiserum by immunoaffinity chromatography. This antibody reacts with human Histone H3. A BLAST analysis was used to suggest cross-reactivity with Human, mouse, and C. elegans. Predicted to react with many species including rat, chicken, Xenopus, Drosophila, and plant based on 100% sequence homology. Cross-reactivity with Histone H3 from other sources has not been determined.
Relevant Links:	• UniProtKB - P84243 • NCBI - NP_002098.1 • GeneID - 3020

Application Details

Tested Applications:	ChIP, Dot Blot, IF, Multiplex, WB
Application Note:	Anti-Histone H3 K79-Me1 antibody is tested in Western Blot, ChIP, Immunofluorescence, and Dot Blot. Specific conditions for reactivity should be optimized by the end user. Expect a band approximately ~15kDa corresponding to the appropriate cell lysate or extract. Epi-Plus™ antibody production in collaboration with Novus Biologicals.
Assay Dilutions:	All assays should be optimized by the user. Recommended dilutions (if any) may be listed below.
ChIP:	2 µg
IF:	1:50 - 1:100
WB:	0.5µg/mL
Other:	Dot Blot - 1:1000

Formulation

Physical State:	Liquid (sterile filtered)
Concentration:	0.50 mg/mL by UV absorbance at 280 nm
Buffer:	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2
Preservative:	0.01% (w/v) Sodium Azide

Stabilizer: 30% Glycerol


Shipping & Handling

Shipping Condition: Dry Ice

Storage Condition: Store vial at -20° C prior to opening. Aliquot contents and freeze at -20° C or below for extended storage. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use.


Expiration: Expiration date is one (1) year from date of receipt.

Images




Immunofluorescence Microscopy

Immunofluorescence of Histone H3 [monomethyl Lys79-Me1]: Histone H3 K79me1 antibody was tested in HeLa cells with DyLight 488 (green). Nuclei and alpha-tubulin were counterstained with DAPI (blue) and DyLight 550 (red).



Dot Blot

Dot Blot of Rabbit Histone H3 [monomethyl Lys79] Antibody. Lane 1: K79 unmodified. Lane 2: K79-Me1. Lane 3: K79-Me2. Lane 4: K79-Me3. Load: 1, 10, and 100 picomoles of peptide. Primary antibody: Histone H3 K79-Me1 antibody at 1:1000 for 45 min at 4°C. Secondary antibody: RABBIT IgG (H&L) Secondary Antibody Peroxidase Conjugated Pre-adsorbed at 1:40,000 for 30 min at RT. Block: 5% BLOTTO 30 minutes at RT.



ChIP

Chromatin Immunoprecipitation of Histone H3 [monomethyl Lys79] Antibody: Chromatin from one million formaldehyde cross-linked HeLa cells was used with 2ug of H3 K79Me1 and 20ul of magnetic IgG beads per immunoprecipitation. A no antibody (No Ab) control was also used. Immunoprecipitated DNA was quantified using quantitative real-time PCR and SYBR green dye, then normalized to the non-precipitated input chromatin, which is equal to one.

Western Blot

Western Blot of Histone H3 [monomethyl Lys79] (RABBIT) Antibody. Western Blot analysis against untreated cell extracts. Lane 1: HeLa cell lysates. Lane 2: NIH/3T3 cell lysates. Lane 3: Cos 7 cell lysates. Load: 35 µg per lane. Primary antibody: Histone H3 K79-Me1 antibody at 0.5µg/mL for overnight at 4°C. Secondary antibody: IRDye800™ rabbit secondary antibody at 1:10,000 for 45 min at RT. Block: 5% BLOTTO overnight at 4°C. Predicted/Observed size: 15 kDa for Histone H3 K79-Me1.

Disclaimer

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 5199, Limerick, Pennsylvania, USA.