

Datasheet for 600-401-P88**Ly6al Antibody****Overview**

Description:	Anti-Ly6al (RABBIT) Antibody - 600-401-P88
Item No.:	600-401-P88
Size:	100 µg
Applications:	WB
Reactivity:	Rat
Host Species:	Rabbit

Product Details

Background:	Ly6al, lymphocyte antigen 6 complex, locus A-like, is also known as Ly6a, Ly6a_predicted. This antibody is suitable for researchers interested in stem cell research.
Synonyms:	Ly6al antibody, lymphocyte antigen 6 complex, locus A-like, Ly6a, Ly6a_predicted, Ly6al antibody
Host Species:	Rabbit
Clonality:	Polyclonal
Format:	IgG

Target Details

Gene Name:	Ly6al
Reactivity:	Rat
Immunogen Type:	Conjugated Peptide
Immunogen:	Ly6al affinity purified antibody was prepared from whole rabbit serum produced by repeated immunizations with a synthetic peptide of rat Ly6al.
Purity/Specificity:	Anti-Lys6al antibody is directed against rat Lys6al protein. The product was affinity purified from monospecific antiserum by immunoaffinity purification. A BLAST analysis was used to suggest reactivity with this protein from rat based on 100% homology for the immunogen sequence. Cross reactivity with Lys6al from other sources has not been determined.

Application Details

Tested Applications:	WB
Application Note:	Anti-Lys6al is tested for Western Blotting. Specific conditions for reactivity should be optimized by the end user. Expect a band approximately ~14.2 kDa corresponding to the appropriate cell lysate or extract.
Assay Dilutions:	All assays should be optimized by the user. Recommended dilutions (if any) may be listed below.
WB:	0.5µg/mL


Formulation

Physical State:	Lyophilized
Concentration:	0.5 mg/mL by UV absorbance at 280 nm
Buffer:	0.9mg NaCl, 0.2mg Na ₂ HPO ₄ , 0.05mg NaN ₃
Preservative:	0.05mg Thimerosal
Stabilizer:	5mg BSA
Reconstitution Volume:	100 µL
Reconstitution Buffer:	Restore with deionized water (or equivalent)

Shipping & Handling

Shipping Condition:	Ambient
Storage Condition:	Store vial at 4° C prior to opening. Aliquot contents and freeze at -20° C or below for extended storage. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use.
Expiration:	Expiration date is one (1) year from date of receipt.

Images

**Western Blot**

Western blotting of Anti-Ly6al antibody.

Lane 1: Rat Cardiac Muscle Tissue Lysate, Lane 2: Rat Liver Tissue Lysate.

Disclaimer

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 5199, Limerick, Pennsylvania, USA.