

Datasheet for 600-401-Z87**B-raf Antibody****Overview**

Description:	Anti-B-raf (RABBIT) Antibody - 600-401-Z87
Item No.:	600-401-Z87
Size:	100 µg
Applications:	ELISA, IF, IHC, WB
Reactivity:	Human, Mouse
Host Species:	Rabbit

Product Details

Background:	B-raf belongs to the raf/mil family of serine/threonine protein kinases and plays a role in regulating the MAP kinase/ERKs signaling pathway, which affects cell division, differentiation, and secretion. The Ras/Raf/MEK/ERK and Ras/PI3K/PTEN/Akt pathways interact with each other to regulate growth and in some cases tumorigenesis. Mutations in B-raf have been associated with several cancers, including non-Hodgkin lymphoma, colorectal cancer, malignant melanoma, thyroid carcinoma, non-small cell lung carcinoma, and adenocarcinoma of lung, leading to speculation on the possibility of B-raf as a therapeutic target for treating cancers. Mutations in this gene have also been associated with cardiofaciocutaneous syndrome (CFCS), a disease characterized by heart defects, mental retardation and a distinctive facial appearance.
Synonyms:	B-raf Antibody, NS7, BRAF1, RAFB1, B-RAF1, Proto-oncogene B-Raf
Host Species:	Rabbit
Clonality:	Polyclonal
Format:	IgG

Target Details

Gene Name:	BRAF
Reactivity:	Human, Mouse
Immunogen Type:	Conjugated Peptide

Immunogen:	Anti-B-raf antibody was prepared from whole rabbit serum produced by repeated immunizations with a 19 amino acid peptide near the internal region of human B-raf. The immunogen is located within amino acids 330 - 380 of B-raf.
Purity/Specificity:	Anti-B-raf Antibody was affinity purified from monospecific antiserum by immunoaffinity chromatography. At least two isoforms of B-raf are known to exist; this antibody will detect both isoforms. This antibody will not cross-react with C-raf. Predicted species reactivity based on immunogen sequence: Rat: (100%), Chicken: (100%).
Relevant Links:	• UniProtKB - P15056 • GeneID - 673 • NCBI - NP_004324.2

Application Details

Tested Applications:	ELISA, IF, IHC, WB
Application Note:	Anti-B-raf Antibody has been tested for use in ELISA, Western Blotting, Immunohistochemistry and Immunofluorescence. Specific conditions for reactivity should be optimized by the end user. Expect a band at approximately 94 kDa in Western Blots of specific cell lysates and tissues. Positive controls used are Human brain tissue, Human Small Intestine Tissue, Daudi cells, and K562 cells.
Assay Dilutions:	All assays should be optimized by the user. Recommended dilutions (if any) may be listed below.
ELISA:	1:10,000-1:20,000
IF:	User Optimized
IHC:	User Optimized
WB:	1 - 2 µg/mL


Formulation

Physical State:	Liquid (sterile filtered)
Concentration:	1 mg/mL by UV absorbance at 280 nm
Buffer:	0.01 M Sodium Phosphate, 0.25 M Sodium Chloride, pH 7.2
Preservative:	0.02% (w/v) Sodium Azide
Stabilizer:	None

Shipping & Handling

Shipping Condition:	Dry Ice
Storage Condition:	Store vial at -20° C prior to opening. Aliquot contents and freeze at -20° C or below for extended storage. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use.
Expiration:	Expiration date is one (1) year from date of receipt.

Images



Immunofluorescence Microscopy


Immunofluorescence of B-raf.

Tissue: Human Small Intestine tissue.

Fixation: 4% paraformaldehyde-fixed.

Primary Antibody: B-raf at 20 µg/mL.

Secondary: goat anti-rabbit IgG secondary antibody at 1:500 dilution (red).



Immunohistochemistry

Immunohistochemistry of B-raf.


Tissue: Human Small Intestine Tissue.

Fixation: paraffin-embedded, fixed with formaldehyde and blocked with 10% serum for 1 h at RT.

Antigen retrieval: heat mediation with a citrate buffer (pH6).

Primary Antibody: anti-B-raf antibody at 2.5 µg/ml overnight at 4°C.

Secondary: goat anti-rabbit IgG H&L (HRP) at 1:250. Counter stained with Hematoxylin.




Western Blot

Western Blot of B-raf.

Load: 15 µg of lysates per lane.

Primary Antibodies: Top: chicken B-raf at 0.5 µg/mL, Middle: rabbit B-raf (p/n 600-401-Z87) at 0.5 µg/mL, and Bottom: GAPDH at 0.02 µg/mL for 1h incubation at RT in 5% NFDm/TBST.

Secondary: Goat anti-chicken or goat anti-rabbit IgG HRP conjugate at 1:10000 dilution.



Western Blot

Western Blot of B-raf.

Load: 15 µg of Human Brain tissue lysates per lane.

Antibodies: B-raf at (A: 1 µg/mL; B: 2 µg/mL) for 1h incubation at RT in 5% NFDm/TBST.

Secondary: Goat anti-rabbit IgG HRP conjugate at 1:10000 dilution.

Disclaimer

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 5199, Limerick, Pennsylvania, USA.