

Datasheet for 600-443-381**MYC Epitope Tag Antibody DyLight™ 649 Conjugated****Overview**

Description:	Anti-Myc Epitope Tag (RABBIT) Antibody DyLight™ 649 Conjugated - 600-443-381
Item No.:	600-443-381
Size:	100 µg
Applications:	WB, Biochemical Assay
Reactivity:	Myc-Tag
Host Species:	Rabbit

Product Details

Background:	Epitope tags are short peptide sequences that are easily recognized by tag-specific antibodies. Due to their small size, epitope tags do not affect the biochemical properties of the tagged protein. The anti-epitope tag antibody is usually functional in a variety of antibody-dependent experimental procedures. Expression vectors producing epitope tag fusion proteins are available for a variety of host expression systems including bacteria, yeast, insect and mammalian cells. Rockland produces anti-epitope tag antibodies against many common epitope tags including Myc, GST, GFP, 6X His, MBP, FLAG and HA. Anti-MYC has utility to detect the fusion protein of the Myc epitope cloned along with the target gene. As such, anti-Myc/Myc can be used to identify fusion proteins containing the Myc epitope. This Anti-Myc antibody has been conjugated to DyLight™649.
Synonyms:	Rabbit anti-MYC Epitope Tag DyLight 649™ Conjugated Antibody, Rabbit anti-c-myc tag DyLight 649™ conjugation, Glu-Gln-Lys-Leu-Ile-Ser-Glu-Glu-Asp-Leu
Host Species:	Rabbit
Conjugate:	DyLight™ 649
Clonality:	Polyclonal
Format:	IgG
F/P Ratio:	1.8

Target Details

Reactivity:	Myc-Tag
--------------------	---------

Immunogen Type:	Conjugated Peptide
Immunogen:	This antibody was purified from whole rabbit serum prepared by repeated immunizations with MYC epitope tag peptide, E-Q-K-L-I-S-E-E-D-L, conjugated to KLH using maleimide. The sequence corresponds to amino acids 410-419 of human c-Myc.
Purity/Specificity:	This affinity purified antibody is directed against the MYC epitope tag and is useful in determining its presence in over expressed proteins in various assays. The antibody recognizes the MYC epitope tag -tag (Glu-Gln-Lys-Leu-Ile-Ser-Glu-Glu-Asp-Leu) fused to either the amino- or carboxy- termini of targeted proteins in transfected or transformed cells.

Application Details

Tested Applications:	WB
Suggested Applications:	Biochemical Assay (Based on references)
Application Note:	Anti-MYC Epitope tag DyLight 649 Conjugated Antibody has been tested by western blot. This product is designed for immunofluorescence microscopy, fluorescence based plate assays (FLISA) and fluorescent western blotting. This product is also suitable for multiplex analysis, including multicolor imaging, utilizing various commercial platforms. The emission spectra for this DyLight™ conjugate match the principle output wavelengths of most common fluorescence instrumentation.
Assay Dilutions:	All assays should be optimized by the user. Recommended dilutions (if any) may be listed below.
FLISA:	>1:20,000
IF:	>1:5,000
WB:	1:10,000 - 1:25,000

Formulation

Physical State:	Lyophilized
Concentration:	1.0 mg/mL by UV absorbance at 280 nm
Buffer:	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2
Preservative:	0.01% (w/v) Sodium Azide
Stabilizer:	10 mg/mL Bovine Serum Albumin (BSA) - Immunoglobulin and Protease free
Reconstitution Volume:	100 µL
Reconstitution Buffer:	Restore with deionized water (or equivalent)


Shipping & Handling

Shipping Condition: Ambient

Storage Condition: Store vial at 4° C prior to restoration. For extended storage aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use.

Expiration: Expiration date is one (1) year from date of receipt.

Images




Figure

KinExA measurements of Axl Ig1 variants. (a) N-curve analysis of Gas6 binding titrations for wild-type Axl Ig1, MYD1 Ig1, and AxlNb Ig1. N-curve analysis of AxlNb could not be performed as complete binding to 1 nM of Gas6 was not attained even at μ M concentrations of AxlNb. This corresponds to a K_D > 1 μ M which is out of the measurable range of this assay. (b) KinExA direct inject kinetic data for wild-type Axl Ig1 and MYD1 Ig1. (c) KinExA signal tests of Protein S and Gas6 binding to MYD1 Fc coated beads. MYD1 Fc does not bind Protein S as no signal was detected after 700 μ l of a 10 nM Protein S solution was passed over beads (left). In comparison, a robust signal was obtained after 100 μ l of a 5 nM solution of Gas6 was flowed over the same beads. Curve fits are shown as dotted black lines and nominal concentrations of Gas6 are indicated. After the appropriate incubation time, reactions were flowed over MYD1 Fc coated beads and captured free Gas6 was probed using an anti-6xHis Dylight 649 antibody (p/n 200-343-382).

To test binding of Protein S to MYD1, 700 μ l of a 10 nM solution was flowed over the capture beads. Binding of Protein S to the beads was probed using an anti-c-Myc Dylight 649 antibody (p/n 600-443-381).

Supplementary Figure 2. PMID: 25242553




Western Blot

DyLight™ dyes can be used for two-color western blot detection with low background and high signal. Anti-tubulin was detected using a DyLight™ 549 conjugate. Anti-TNFα was detected using a DyLight™ 649 conjugate. The image was captured using the Typhoon™ 9410 Imaging System.

Diagram


Properties of DyLight™ Fluorescent Dyes.

Emission	Color	DyLight™ Dye	Ex/Em (nm)	ϵ (M ⁻¹ cm ⁻¹)	Similar Dyes
Blue		405	400/420	30,000	Alexa™ 405, Cascade Blue
Green		488	493/518	70,000	Alexa™ 488, Cy2®, FITC
Yellow		549	550/568	150,000	Alexa™ 546, Alexa 555, Cy3®, TRITC
Red		649	646/674	250,000	Alexa™ 647, Cy5®
Near Infrared		680	682/715	140,000	Alexa™ 680, Cy5.5®, IRDye™ 700
Infrared		800	770/794	270,000	IRDye™ 800



Western Blot

Western Blot of Rabbit Anti-MYC Antibody DyLight™649 Conjugated. Lane 1: Opal Prestained Molecular Weight Marker (p/n MB-210-0500). Lane 2: Flag Lysate (p/n W00-001-383). Lane 3: 12 Epitope GST Tag Protein Marker Lysate (p/n MB-302-0100). Lane 4: HeLa Whole Cell Lysate (p/n W09-000-364). Load: 10µg/lane. Primary Antibody: Anti-Myc Epitope Tag DyLight™649 at 1.0µg/mL overnight at 2-8°C. Blocking: MB-070 Buffer at RT for 1hr. Expected MW: ~64kDa for Flag Lysate, ~48kDa for 12 Epitope GST Tag. Exposure: 2.2 sec.



Diagram

References

- Kariolis MS et al. An engineered Axl decoy receptor effectively silences the Gas6-Axl signaling axis. *Nat Chem Biol.* (2014)

Disclaimer

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 5199, Limerick, Pennsylvania, USA.