

Datasheet for 605-309-002

## Goat IgG (H&L) Antibody Texas Red Conjugated

### Overview

|                      |   |
|----------------------|---|
| <b>Description:</b>  | Mouse Anti-Goat IgG (H&L) Antibody Texas Red Conjugated - 605-309-002 |
| <b>Item No.:</b>     | 605-309-002   |
| <b>Size:</b>         | 1 mg  |
| <b>Applications:</b> | Dot Blot  |
| <b>Reactivity:</b>   | Goat  |
| <b>Host Species:</b> | Mouse   |

### Product Details

|                      |   |
|----------------------|---|
| <b>Background:</b>   | Secondary Antibodies are available in a variety of formats and conjugate types. When choosing a secondary antibody product, consideration must be given to species and immunoglobulin specificity, conjugate type, fragment and chain specificity, level of cross-reactivity, and host-species source and fragment composition. |
| <b>Synonyms:</b>     | Mouse anti-Goat IgG Antibody Texas Red™ Conjugation, Mouse anti-Goat IgG Texas Red-™ Conjugated Antibody  |
| <b>Host Species:</b> | Mouse   |
| <b>Specificity:</b>  | IgG (H&L)   |
| <b>Conjugate:</b>    | Texas Red®  |
| <b>Clonality:</b>    | Polyclonal  |
| <b>Format:</b>       | IgG   |
| <b>F/P Ratio:</b>    | 2.3   |

### Target Details

|                    |                         |
|--------------------|-------------------------|
| <b>Reactivity:</b> | Goat                    |
| <b>Immunogen:</b>  | Goat IgG whole molecule |

**Purity/Specificity:** Goat IgG (H&L) Antibody Texas Red Conjugated was prepared from monospecific polyclonal ascites by immunoaffinity chromatography using Goat IgG coupled to agarose beads followed by solid phase adsorption(s) to remove any unwanted reactivities. Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-Texas Red, anti-Mouse Serum, Goat IgG and Goat Serum. No reaction was observed against Rabbit, Human and Mouse Serum Proteins.

## Application Details

**Tested Applications:** Dot Blot

**Application Note:** Mouse Anti-Goat IgG Texas Red Conjugate has been tested by dot blot and is designed for immunofluorescence microscopy, fluorescence based plate assays (FLISA) and fluorescent western blotting. This product is also suitable for multiplex analysis, including multicolor imaging, utilizing various commercial platforms.

**Assay Dilutions:** All assays should be optimized by the user. Recommended dilutions (if any) may be listed below.

**IF:** 1:500-1:2,500

## Formulation

**Physical State:** Lyophilized

**Concentration:** 1.0 mg/mL by UV absorbance at 280 nm

**Buffer:** 0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2

**Preservative:** 0.01% (w/v) Sodium Azide

**Stabilizer:** 10 mg/mL Bovine Serum Albumin (BSA) - Immunoglobulin and Protease free

**Reconstitution Volume:** 1.0 mL

**Reconstitution Buffer:** Restore with deionized water (or equivalent)

## Shipping & Handling

**Shipping Condition:** Ambient

**Storage Condition:** Store vial at 4° C prior to restoration. For extended storage aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use.

**Expiration:** Expiration date is one (1) year from date of receipt.

## Disclaimer

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 5199, Limerick, Pennsylvania, USA.