

## Datasheet for 611-104-122

**Rabbit IgG (H&L) Antibody CY3 Conjugated Pre-adsorbed****Overview**

<b>Description:</b>	Goat Anti-Rabbit IgG (H&L) Antibody CY3 Conjugated (Min X Bv Ch Gt GP Ham Hs Hu Ms Rt & Sh Serum Proteins) - 611-104-122
<b>Item No.:</b>	611-104-122
<b>Size:</b>	1 mg
<b>Applications:</b>	Dot Blot, WB, IF, IHC, Multiplex
<b>Reactivity:</b>	Rabbit
<b>Host Species:</b>	Goat

**Product Details**

**Background:** Anti-Rabbit IgG CY3 Conjugated Antibody generated in goat detects rabbit IgG. Secreted as part of the adaptive immune response by plasma B cells, immunoglobulin G constitutes 75% of serum immunoglobulins. Immunoglobulin G binds to viruses, bacteria, as well as fungi and facilitates their destruction or neutralization via agglutination (and thereby immobilizing them), activation of the complement cascade, and opsonization for phagocytosis. The whole IgG molecule possesses both the F(c) region, recognized by high-affinity Fc receptor proteins, as well as the F(ab) region possessing the epitope-recognition site. Both heavy and light chains of the antibody molecule are present. Secondary Antibodies are available in a variety of formats and conjugate types. When choosing a secondary antibody product, consideration must be given to species and immunoglobulin specificity, conjugate type, fragment and chain specificity, level of cross-reactivity, and host-species source and fragment composition. This Anti-Rabbit IgG (H&L) is conjugated to Cy3.

<b>Synonyms:</b>	goat anti-Rabbit IgG CY3 Conjugated antibody, goat anti-Rabbit IgG antibody CY3 conjugation
<b>Host Species:</b>	Goat
<b>Specificity:</b>	IgG (H&L)
<b>Conjugate:</b>	Cy3™
<b>Clonality:</b>	Polyclonal
<b>Format:</b>	IgG
<b>F/P Ratio:</b>	10.1

## Target Details

<b>Reactivity:</b>	Rabbit
<b>Immunogen:</b>	Rabbit IgG whole molecule
<b>Purity/Specificity:</b>	RABBIT IgG (H&L) Antibody CY3 Conjugated Pre-adsorbed was prepared from monospecific antiserum by immunoaffinity chromatography using Rabbit IgG coupled to agarose beads followed by solid phase adsorption(s) to remove any unwanted reactivities and extensive dialysis against the buffer stated above. Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-Goat Serum, Rabbit IgG and Rabbit Serum. No reaction was observed against Bovine, Chicken, Goat, Guinea Pig, Hamster, Horse, Human, Mouse, Rat and Sheep Serum Proteins.

## Application Details

<b>Tested Applications:</b>	Dot Blot, WB
<b>Suggested Applications:</b>	IF, IHC, Multiplex (Based on references)
<b>Application Note:</b>	Anti-Rabbit IgG (H&L) Antibody CY3 Conjugated Pre-adsorbed antibodies has been tested by dot blot and western blot and are designed for immunofluorescence microscopy, fluorescence based plate assays (FLISA) and fluorescent western blotting. This product is also suitable for multiplex analysis, including multicolor imaging, utilizing various commercial platforms.
<b>Assay Dilutions:</b>	All assays should be optimized by the user. Recommended dilutions (if any) may be listed below.
<b>FC:</b>	1:500 - 1:2,500
<b>FLISA:</b>	1:10,000 - 1:50,000
<b>IF:</b>	1:1,000 - 1:5,000
<b>WB:</b>	User Optimized

## Formulation

<b>Physical State:</b>	Lyophilized
<b>Concentration:</b>	1.0 mg/mL by UV absorbance at 280 nm
<b>Buffer:</b>	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2
<b>Preservative:</b>	0.01% (w/v) Sodium Azide
<b>Stabilizer:</b>	10 mg/mL Bovine Serum Albumin (BSA) - Immunoglobulin and Protease free

**Reconstitution Volume:** 1.0 mL

**Reconstitution Buffer:** Restore with deionized water (or equivalent)


## Shipping & Handling

**Shipping Condition:** Ambient

**Storage Condition:** Store vial at 4° C prior to restoration. For extended storage aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use.


**Expiration:** Expiration date is one (1) year from date of receipt.

## Images




### Dot Blot

Dot Blot Results of Goat Anti-Rabbit IgG CY3 Conjugated Antibody Min X Bv Ch Gt GP Ham Hs Hu Ms Rt & Sh Serum Proteins. Rabbit IgG (p/n 011-0102) samples 1-100ng, 2-33.3ng, 3- 11.1ng, 4- 3.70ng, 5- 1.23ng. Goat Anti-Rabbit IgG CY3 Conjugated Antibody at 1.0µg/mL for 30mins at RT. Blocking Buffer for Fluorescent Western Blotting (p/n MB-070) for 60mins at RT.



#### Western Blot

Western Blot of Unconjugated Anti-Rabbit IgG (H&L) (GOAT) Antibody (Min X Bv, Ch, Gt, GP, Ham, Hs, Hu, Ms, Rt & Sh Serum Proteins) (p/n 611-101-122). Lane M: 3  $\mu$ l Molecular Ladder. Lane 1: Rabbit IgG whole molecule (p/n 011-0102). Lane 2: Rabbit IgG F(ab) Fragment (p/n 011-0105). Lane 3: Rabbit IgG F(c) Fragment (p/n 010-0103). Lane 4: Rabbit IgM Whole Molecule (p/n 011-0107). Lane 5: Normal Rabbit Serum (p/n B309). All samples were reduced. Load: 50 ng per lane. Block: MB-070 for 30 min at RT. Primary Antibody: Anti-Rabbit IgG (H&L) (GOAT) Antibody (Min X Bv, Ch, Gt, GP, Ham, Hs, Hu, Ms, Rt & Sh Serum Proteins) (p/n 611-101-122) 1:1,000 for 60 min at RT. Secondary antibody: Anti-Goat IgG (DONKEY) Peroxidase Conjugated Antibody (p/n CUST10) 1:40,000 in MB-070 for 30 min at RT. Predicted/Observed Size: 25 and 50 kDa for Rabbit IgG and Serum, 25 kDa for F(c) and F(ab), 70 and 23 kDa for IgM. Rabbit F(c) migrates slightly higher.



#### Western Blot

Western Blot of Goat anti-Rabbit IgG Pre-Absorbed Cy3 Conjugated Secondary Antibody. Lane 1: Rabbit IgG. Load: 50 ng per lane. Primary antibody: None. Secondary antibody: Cy3 goat anti-rabbit secondary antibody at 1:1,000 for 60 min at RT. Block: MB-070 for 30 min at RT. Predicted/Observed size: 28 & 55 kDa, 28 & 55 kDa for Rabbit IgG. Other band(s): None.

## References

- Varlamova EG et al. Generation and Characterization of Three Novel Mouse Mutant Strains Susceptible to Audiogenic Seizures. *Cells*. (2024)
- Tao, X et al. Phenylbutyrate modulates polyamine acetylase and ameliorates Snyder-Robinson syndrome in a Drosophila model and patient cells. *Jci Insight* (2022)
- Rubel D et al. Anti-microRNA-21 Therapy on Top of ACE Inhibition Delays Renal Failure in Alport Syndrome Mouse Models. *Cells*. (2022)
- Ma X et al. Nicotinamide mononucleotide adenylyltransferase uses its NAD<sup>+</sup> substrate-binding site to chaperone phosphorylated Tau. *Elife*. (2020)
- Acosta AJ et al. Phosphorylation of GATA4 at serine 105 is required for left ventricular remodelling process in angiotensin II-induced hypertension in rats. *Basic Clin Pharmacol Toxicol*. (2020)
- Rubel D et al. Antifibrotic, nephroprotective effects of paricalcitol versus calcitriol on top of ACE-inhibitor therapy in the COL4A3 knockout mouse model for progressive renal fibrosis. *Nephrol Dial Transplant*. (2014)
- Hartmann K et al. Amniotic fluid derived stem cells give rise to neuron-like cells without a further differentiation potential into retina-like cells. *Am J Stem Cells*. (2013)
- Yin JY et al. Study on the method of quantitative analysis of serum ferritin and soluble transferrin receptor with protein microarray technology. *Biomed Environ Sci*. (2012)
- Ou Q et al. Nuclear receptor DHR4 controls the timing of steroid hormone pulses during Drosophila development. *PLoS One*. (2011)
- Arnhold S et al. Neuronal characteristics of amniotic fluid derived cells after adenoviral transformation. *Cell Biol Int*. (2008)
- Outeiro TF et al. Formation of toxic oligomeric alpha-synuclein species in living cells. *PLoS One* (2008)
- Schurmann FW et al. Glutamate-like immunoreactivity marks compartments of the mushroom bodies in the brain of the cricket. *J Comp Neurol*. (2000)
- Wolter R et al. Mapping of cytoskeletal components in the hippocampal formation of the tree shrew (*Tupaia belangeri*). *J Chem Neuroanat*. (1999)

## Disclaimer

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 5199, Limerick, Pennsylvania, USA.