

## Datasheet for 614-4102

**Swine IgG (H&L) Antibody****Overview**

|                      |                                                 |
|----------------------|-------------------------------------------------|
| <b>Description:</b>  | Rabbit Anti-Swine IgG (H&L) Antibody - 614-4102 |
| <b>Item No.:</b>     | 614-4102                                        |
| <b>Size:</b>         | 10 mg                                           |
| <b>Applications:</b> | ELISA                                           |
| <b>Reactivity:</b>   | Pig                                             |
| <b>Host Species:</b> | Rabbit                                          |

**Product Details**

|                      |                                                                                                                                                                                                       |
|----------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Background:</b>   | SWINE IgG (H&L) (RABBIT) Antibody generated in rabbit detects specifically swine IgG whole molecule. Anti-Swine IgG antibody is ideal for investigators involved in serum protein component research. |
| <b>Synonyms:</b>     | Rabbit anti-Swine Antibody, Anti-Pig IgG (H&L) (Rabbit), Swine Secondary Antibody, Swine in Rabbit Ab, Rabbit anti Swine Antibody                                                                     |
| <b>Host Species:</b> | Rabbit                                                                                                                                                                                                |
| <b>Specificity:</b>  | IgG (H&L)                                                                                                                                                                                             |
| <b>Clonality:</b>    | Polyclonal                                                                                                                                                                                            |
| <b>Format:</b>       | IgG                                                                                                                                                                                                   |

**Target Details**

|                            |                                                                                                                                                                                                                                                                                                                                                                |
|----------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Reactivity:</b>         | Pig                                                                                                                                                                                                                                                                                                                                                            |
| <b>Immunogen:</b>          | Swine IgG whole molecule                                                                                                                                                                                                                                                                                                                                       |
| <b>Purity/Specificity:</b> | Anti-SWINE IgG (H&L) (RABBIT) Antibody was prepared from monospecific antiserum by immunoaffinity chromatography using Swine IgG coupled to agarose beads followed by solid phase adsorption(s) to remove any unwanted reactivities. Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-Rabbit Serum, Swine IgG, and Swine Serum. |

## Application Details

|                             |                                                                                                                                                                                                                                                                                                                                               |
|-----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Tested Applications:</b> | ELISA                                                                                                                                                                                                                                                                                                                                         |
| <b>Application Note:</b>    | Anti-SWINE IgG (H&L) (RABBIT) Antibody has been tested by ELISA and is suitable for immunoblotting (western or dot blot), ELISA, immunoperoxidase electron microscopy and immunohistochemistry as well as other peroxidase-antibody based enzymatic assays requiring lot-to-lot consistency. Specific conditions should be optimized by user. |
| <b>Assay Dilutions:</b>     | All assays should be optimized by the user. Recommended dilutions (if any) may be listed below.                                                                                                                                                                                                                                               |
| <b>ELISA:</b>               | 1:20,000 - 1:100,000                                                                                                                                                                                                                                                                                                                          |
| <b>IHC:</b>                 | 1:1,000 - 1:5,000                                                                                                                                                                                                                                                                                                                             |
| <b>WB:</b>                  | 1:2,000 - 1:10,000                                                                                                                                                                                                                                                                                                                            |


## Formulation

|                        |                                                         |
|------------------------|---------------------------------------------------------|
| <b>Physical State:</b> | Liquid (sterile filtered)                               |
| <b>Concentration:</b>  | 8.8mg/mL by UV absorbance at 280 nm                     |
| <b>Buffer:</b>         | 0.01 M Sodium Phosphate, 0.15 M Sodium Chloride, pH 7.2 |
| <b>Preservative:</b>   | 0.01% (w/v) Sodium Azide                                |
| <b>Stabilizer:</b>     | None                                                    |

## Shipping & Handling

|                            |                                                                                                                                                                                                                                                              |
|----------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Shipping Condition:</b> | Wet Ice                                                                                                                                                                                                                                                      |
| <b>Storage Condition:</b>  | Store vial at 4° C prior to opening. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. |
| <b>Expiration:</b>         | Expiration date is one (1) year from date of receipt.                                                                                                                                                                                                        |

## Images



### ELISA

ELISA Results of Rabbit Anti-Swine IgG Antibody tested against purified Swine IgG. Each well was coated in duplicate with 1.0 µg of Swine IgG (p/n 014-0102). The working dilution is 1:607,000. The starting dilution of antibody was 5µg/ml and the X-axis represents the Log10 of a 3-fold dilution. This titration is a 4-parameter curve fit where the IC50 is defined as the titer of the antibody. Assay performed using HRP Conjugation Stabilizer (p/n MB-076), Goat Anti-Rabbit IgG HRP conjugated (p/n 611-103-122), and TMB substrate (p/n TMBE-1000).

## References

- To H et al. Application of an enzyme-linked immunosorbent assay to detect antibodies to *Actinobacillus pleuropneumoniae* serovars 1, 2 and 5 in the pig sera. *J Vet Med Sci.* (2026)

## Disclaimer

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 5199, Limerick, Pennsylvania, USA.