

Datasheet for 700-109-096

## F(Ab')<sub>2</sub> Fluorescein Antibody Texas Red™ Conjugated

### Overview

<b>Description:</b>	F(ab') <sub>2</sub> Anti-Fluorescein (GOAT) Antibody Texas Red™ Conjugated - 700-109-096
<b>Item No.:</b>	700-109-096
<b>Size:</b>	1 mg
<b>Reactivity:</b>	Fluorescein
<b>Host Species:</b>	Goat

### Product Details

<b>Background:</b>	F(ab') <sub>2</sub> Anti-Fluorescein Antibody generated in goat detects fluorescein. It is an organic dye used as a fluorescent tracer in many application including cellular staining. This F(ab') <sub>2</sub> Anti-Fluorescein Antibody is Texas Red™ conjugated.
<b>Synonyms:</b>	goat F(ab') <sub>2</sub> anti-Fluorescein Antibody Texas Red™ Conjugation, goat F(ab') <sub>2</sub> anti-Fluorescein Texas Red™ Conjugated Antibody, goat F(ab') <sub>2</sub> anti-FITC Antibody TR Conjugation
<b>Host Species:</b>	Goat
<b>Conjugate:</b>	Texas Red®
<b>Clonality:</b>	Polyclonal
<b>Format:</b>	IgG F(ab') <sub>2</sub>
<b>F/P Ratio:</b>	3.5

### Target Details

<b>Reactivity:</b>	Fluorescein
<b>Immunogen Type:</b>	Other
<b>Immunogen:</b>	Fluorescein conjugated to Goat IgG
<b>Purity/Specificity:</b>	This product was prepared from monospecific antiserum by immunoaffinity chromatography using Fluorescein coupled to agarose beads followed by solid phase adsorption(s) to remove any unwanted reactivities, pepsin digestion and chromatographic separation. Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-Goat Serum and Fluorescein conjugated BSA. No reaction was observed against anti-Pepsin or anti-Goat IgG F(c).

## Application Details

<b>Application Note:</b>	This product is designed for immunofluorescence microscopy, fluorescence based plate assays (FLISA) and fluorescent western blotting. This product is also suitable for multiplex analysis, including multicolor imaging, utilizing various commercial platforms as well as other antibody based fluorescent assays requiring lot-to-lot consistency.
<b>Assay Dilutions:</b>	All assays should be optimized by the user. Recommended dilutions (if any) may be listed below.
<b>FLISA:</b>	>1:20,000
<b>IF:</b>	1:500 - 1:2,500
<b>WB:</b>	>1:10,000

## Formulation

<b>Physical State:</b>	Lyophilized
<b>Concentration:</b>	1.0 mg/mL by UV absorbance at 280 nm
<b>Buffer:</b>	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2
<b>Preservative:</b>	0.01% (w/v) Sodium Azide
<b>Stabilizer:</b>	10 mg/mL Bovine Serum Albumin (BSA) - Immunoglobulin and Protease free
<b>Reconstitution Volume:</b>	1.0 mL
<b>Reconstitution Buffer:</b>	Restore with deionized water (or equivalent)

## Shipping & Handling

<b>Shipping Condition:</b>	Ambient
<b>Storage Condition:</b>	Store vial at 4° C prior to restoration. For extended storage aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use.
<b>Expiration:</b>	Expiration date is one (1) year from date of receipt.

## Disclaimer

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 5199, Limerick, Pennsylvania, USA.