

Datasheet for B309**Normal Rabbit Serum (NRS)****Overview**

Description:	Normal Rabbit Serum (NRS) - B309
Item No.:	B309
Size:	10 mL
Applications:	IHC, WB
Origin:	Rabbit

Product Details

Background:	Normal Rabbit Serum (NRS) is the component of blood that does not contain erythrocytes or leukocytes nor clotting factors. All other non-clotting proteins remain in the serum making it a useful reagent for blocking of blots or wells. Normal Rabbit Serum is non-sterile and not amenable to cell culture.
Synonyms:	normal rabbit serum, blocking buffer, NRS, Rabbit sera, non-sterile rabbit sera, normal sera
Species of Origin:	Rabbit

Target Details

Purity/Specificity:	Normal Rabbit Serum (NRS) was prepared from normal serum by a multi-step process which includes delipidation and selective precipitation. Normal Rabbit Serum (NRS) assayed by immunoelectrophoresis resulted in a multiple precipitin arcs against anti-Rabbit Serum.
----------------------------	--

Application Details

Suggested Applications:	IHC, WB (Based on references)
Application Note:	Normal Rabbit Serum (NRS) is a useful reagent for blocking Western Blots and ELISA wells.
Assay Dilutions:	All assays should be optimized by the user. Recommended dilutions (if any) may be listed below.
ELISA:	1x
WB:	1x

Tissue Data

Tissue Type:	Serum
---------------------	-------


Formulation

Physical State:	Lyophilized
Concentration:	90 mg/mL by Refractometry
Buffer:	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2
Sterility:	Non-sterile
Preservative:	0.01% (w/v) Sodium Azide
Reconstitution Volume:	10 mL
Reconstitution Buffer:	Restore with deionized water (or equivalent)

Shipping & Handling

Shipping Condition:	Ambient
Storage Condition:	Store vial at 4° C prior to restoration. For extended storage aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. Normal Rabbit Serum (NRS) is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use.
Expiration:	Expiration date is one (1) year from date of receipt.

Images



Bottle
Normal Rabbit Serum (NRS)

References

- Cuoghi I et al. Immunohistochemical analysis of the distribution of molecules involved in ionic and pH regulation in the lancelet *Branchiostoma floridae* (Hubbs, 1922). *Acta Histochem.* (2018)
- Mohawk, KL et al. Neutralizing antibodies to Shiga toxin type 2 (Stx2) reduce colonization of mice by Stx2-expressing *Escherichia coli* O157:H7. *Vaccine* (2010)

Disclaimer

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 5199, Limerick, Pennsylvania, USA.