

Datasheet for C201-0005**Mouse Complement (Lyophilized) With Diluent****Overview**

Description:	Mouse Complement (Lyophilized) With Diluent - C201-0005
Item No.:	C201-0005
Size:	5 mL
Applications:	Other
Origin:	Mouse

Product Details

Background: The complement system helps or “complements” the ability of antibodies and phagocytic cells to clear pathogens from an organism. Complement is part of the innate immune system that is not adaptable and does not change over the course of an individual's lifetime. Complement can be recruited and brought into action by the adaptive immune system.

The following are basic functions of complement:

1. Opsonization - enhancing phagocytosis of antigens
2. Chemotaxis - attracting macrophages and neutrophils
3. Cell Lysis - rupturing membranes of foreign cells
4. Clumping of antigen-bearing agents

Synonyms: Complement system, tissue macrophages, blood monocytes, protease C3-convertase, mannose-binding lectin pathway, C3, C3a, C3b, C5a, C5b, C6, C7, C8, and polymeric C9, cascade cleavage and activation events, recruit inflammatory cells, anaphylatoxin

Species of Origin: Mouse

Target Details

Purity/Specificity: The mouse complement system consists of a number of small proteins found in the blood, generally synthesized by the liver, and normally circulating as inactive precursors (pro-proteins). When stimulated by one of several triggers, proteases in the system cleave specific proteins to release cytokines and initiate an amplifying cascade of further cleavages. The end-result of this activation cascade is massive amplification of the response and activation of the cell-killing membrane attack complex. Over 25 proteins and protein fragments make up the complement system, including serum proteins, serosal proteins, and cell membrane receptors. They account

for about 5% of the globulin fraction of blood serum.

Application Details

Suggested Applications:	Other (Based on references)
Application Note:	pH: normal Immunoelectrophoresis: normal Hemoglobin: normal IgG Concentration: normal Immunoelectrophoresis: normal
Assay Dilutions:	All assays should be optimized by the user. Recommended dilutions (if any) may be listed below.

Tissue Data

Tissue Type:	Complement
Sex:	Mixed
Strain:	Mouse - Mixed

Formulation

Physical State:	Lyophilized
Concentration:	84 mg/ml by Refractometry
Buffer:	None
Sterility:	Non-sterile
Preservative:	None
Stabilizer:	None
Reconstitution Volume:	5.0 mL
Reconstitution Buffer:	Restore with Diluent

Shipping & Handling

Shipping Condition:	Ambient
Storage Condition:	Store vial at 4° C prior to restoration. Aliquot contents and freeze at -70° C or below. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. COMPLEMENT IS A TEMPERATURE SENSITIVE PRODUCT. IMPROPER STORAGE WILL INACTIVATE COMPLEMENT ACTIVITY.
Expiration:	Expiration date is one (1) year from date of receipt.

References

- Farías MA et al. Preclinical characterisation of the protective capacity of an anti-nucleoprotein hRSV monoclonal antibody. *EBioMedicine*. (2026)
- Hijaze N et al. Inducing regulated necrosis and shifting macrophage polarization with anti-EMMPRIN antibody (161-pAb) and complement factors. *J Leukoc Biol*. (2020)

Disclaimer

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 5199, Limerick, Pennsylvania, USA.