

Datasheet for C220-0002

Rabbit Anti-Sheep Cell Hemolysin

Overview

Description:	Sheep Cell Hemolysin (RABBIT) Antibody - C220-0002
Item No.:	C220-0002
Size:	2 mL
Applications:	Agglutination, ELISA
Reactivity:	Sheep
Host Species:	Rabbit

Product Details

Background:	Sheep red blood cells coated with Rabbit anti-Sheep Cell Hemolysin have traditionally been used to assay the activity of the complement classical pathway in serum samples.
Synonyms:	Rabbit anti-sheep erythrocyte antiserum, Anti-RBC antibody, Red Blood Cell Antibody, Antibody for hemagglutination, rabbit anti-sheep cell hemolysin, rabbit anti-sheep RBC, rabbit antibody to sheep Red Blood Cells (RBC), haemolysin, hemolysin, erythrocytes sensitizing agent, sheep hemolysin
Host Species:	Rabbit
Clonality:	Polyclonal

Target Details

Reactivity:	Sheep
Immunogen Type:	Other
Immunogen:	Sheep washed pooled Red Blood Cells (RBC)
Purity/Specificity:	Rabbit Sheep Cell Hemolysin was prepared from polyspecific antiserum by a delipidation and defibrination.
Relevant Links:	C220 SDS

Application Details

Tested Applications:	Agglutination
Suggested Applications:	ELISA (Based on references)
Application Note:	Rabbit Anti-Sheep Cell Hemolysin has been tested by agglutination and is suitable for agglutination, complement fixation related assays.
Assay Dilutions:	All assays should be optimized by the user. Recommended dilutions (if any) may be listed below.

Tissue Data

Tissue Type:	Red Blood Cells
---------------------	-----------------


Formulation

Physical State:	Lyophilized
Concentration:	75 mg/mL by Refractometry
Buffer:	None
Sterility:	Non-sterile
Preservative:	None
Stabilizer:	None
Reconstitution Volume:	2.0 mL
Reconstitution Buffer:	Restore with deionized water (or equivalent)

Shipping & Handling

Shipping Condition:	Ambient
Storage Condition:	Store vial at 4° C prior to restoration. For extended storage aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use.
Expiration:	Expiration date is one (1) year from date of receipt.

Images



ELISA

Haemolytic activity of disulphide locked C9 against erythrocytes/antibody/complement 1–8 (EAC1-8). The TMH1 locked (no DTT) alone is inactive; however, activity can be rescued with 1 mM DTT (TMH1 locked (with DTT)). Also shown are control experiments: no C9, and wildtype C9 (with and without DTT). EAC1–8 was produced by sensitising 6.5×10^8 cells mL⁻¹ sRBC with equal volume of anti-sheep antibody (0.75 mg mL⁻¹) (Rockland Immunochemicals p/n C220–0002) at 30 °C. Fig 2. PMID: 30111885.

References

- Spicer BA et al. The first transmembrane region of complement component-9 acts as a brake on its self-assembly. *Nat Commun.* (2018)

Disclaimer

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 5199, Limerick, Pennsylvania, USA.