

**Datasheet for DAB-B50****10X Concentrate Diluent Buffer for DAB Membrane Substrate****Overview**

|                     |                                                                              |
|---------------------|------------------------------------------------------------------------------|
| <b>Description:</b> | 10X Concentrate Diluent Buffer For DAB Membrane Substrate (DAB-10) - DAB-B50 |
| <b>Item No.:</b>    | DAB-B50                                                                      |
| <b>Size:</b>        | 50 mL                                                                        |

**Product Details**

|                  |                                                           |
|------------------|-----------------------------------------------------------|
| <b>Synonyms:</b> | Dab Membrane Peroxidase Diluent, Dab Membrane HRP Diluent |
|------------------|-----------------------------------------------------------|

**Target Details**

|                            |                                                 |
|----------------------------|-------------------------------------------------|
| <b>Purity/Specificity:</b> | pH: 6.0 +/- 0.1                                 |
|                            | Stability @ 37°C: PASS                          |
|                            | Stability @ 18° to 26°C: PASS                   |
|                            | Stability @ 4°C: PASS                           |
|                            | QC Raw Material: PASS                           |
|                            | Absorbance check of final product: PASS         |
|                            | Performance Data per Immunohistochemistry: PASS |

**Application Details**

|                          |                                                                                                 |
|--------------------------|-------------------------------------------------------------------------------------------------|
| <b>Application Note:</b> | Recommended to be used with DAB-10.                                                             |
| <b>Assay Dilutions:</b>  | All assays should be optimized by the user. Recommended dilutions (if any) may be listed below. |
| <b>IHC:</b>              | 1X                                                                                              |
| <b>WB:</b>               | 1X                                                                                              |

## Formulation

|                        |                                                    |
|------------------------|----------------------------------------------------|
| <b>Physical State:</b> | Liquid - clear                                     |
| <b>Stabilizer:</b>     | Proprietary 10X TRIS buffer with hydrogen peroxide |

## Shipping & Handling

|                            |                                                                                                                                        |
|----------------------------|----------------------------------------------------------------------------------------------------------------------------------------|
| <b>Shipping Condition:</b> | Ambient                                                                                                                                |
| <b>Storage Condition:</b>  | Store container at 4° C prior to opening. Protect from moisture and light. No special shipping conditions or precautions are required. |
| <b>Expiration:</b>         | Expiration date is one (1) year from date of receipt.                                                                                  |

## Images



### Bottle

10X Concentrate Diluent Buffer For DAB Membrane Substrate (DAB-10)

## Disclaimer

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 5199, Limerick, Pennsylvania, USA.

