

Datasheet for GP-T038

Guinea Pig Adrenal

Overview

Description:	Guinea Pig Adrenal - GP-T038
Item No.:	GP-T038
Size:	1 Each
Origin:	Guinea Pig

Product Details

Species of Origin:	Guinea Pig
--------------------	------------

Formulation

Physical State:	Tissue
Sterility:	Non-sterile

Shipping & Handling

Shipping Condition:	Dry Ice
Storage Condition:	Store tissue at -20° C or colder prior to use.
Expiration:	No expiration date is given for this product if properly stored.

Disclaimer

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 5199, Limerick, Pennsylvania, USA.

