

Datasheet for KCAJ07**E. coli HCP Chemiluminescent Western Blot Kit****Overview**

Description:	Chemiluminescent Western Blot Kit for detection of E.coli Host Cell Proteins (HCP) - KCAJ07
Item No.:	KCAJ07
Size:	1 Kit
Applications:	WB
Reactivity:	E. coli

Product Details

Background:	Host Cell Proteins are process-related impurities derived from a host cell expression system that may be present in trace amounts in a final drug substance. HCP Antibodies are used to detect impurities in biologically-based therapeutics. The antibody provided in this kit recognizes high and low molecular weight protein fractions from E.coli. The Anti-Host Cell Proteins Antibody provided in this kit detects over 300 E.coli proteins and is suitable for western blotting. This antibody has been validated by 2D electrophoresis western blot to assure highest possible coverage of E.coli HCP. Detecting impurities often unseen by 1-D electrophoresis, a 2D electrophoresis western blot analysis delivers more complete insight into the immunocoverage of protein components of the host cell lysate.
Synonyms:	Western Blotting Kit, Chemiluminescent Kit, Peroxidase Kit, Immunoblotting Kit, HCP, Host Cell Protein detection, HCP contaminants detection, ECP, ECP detection, ECP contaminants
Detection Kit Type:	Chemiluminescent Western Blot Kit

Target Details

Reactivity:	E. coli
Immunogen:	Anti-Host Cell Protein Antibody (p/n 100-401-J07) was prepared blending antibodies that have been separately developed by repeated immunizations in rabbit with a mixture of high molecular weight (HMW) and a mixture of low molecular weight (LMW) proteins from E.Coli (BL21).
Purity/Specificity:	This Chemiluminescent Western Blot Kit allows for the detection of E.coli host cell proteins present in samples provided by the user. Included in the kit along with Rabbit Anti-Host Cell Protein antibody is a positive protein control, Goat Anti-Rabbit IgG HRP Antibody, and Blocking Buffer.

Relevant Links:

- [Chemiluminescent Western Blot Kit for detection of E.coli HCP - KCAJ07](#)

Application Details

Tested Applications:	WB
Application Note:	Use Rockland Immunochemicals' Chemiluminescent Western Blot Kit for detection of E.coli Host Cell Proteins. This kit is useful for both "Western blotting" and "dot blotting" methods. After protein separation and transfer, the membrane is probed with Anti-Host Cell Proteins Antibody. Detection of the membrane bound antibody-antigen complex is achieved by the addition of a secondary antibody conjugated to the enzyme horseradish peroxidase. The enzyme reacts with Femtomax, a specialized formulation of luminol which is an extremely sensitive, non-radioactive substrate that emits light (chemiluminescence) and allows visualization using X-ray film or other imaging methods, including highly sensitive CCD cameras and imaging systems. Femtomax reagent (FEMTOMAX-020 or FEMTOMAX-110) is not included with this kit and must be purchased separately.
Assay Dilutions:	All assays should be optimized by the user. Recommended dilutions (if any) may be listed below.
WB:	User Optimized by Kit Protocol


Formulation

Physical State:	n/a
------------------------	-----

Shipping & Handling

Shipping Condition:	Wet Ice
Storage Condition:	See kit insert for complete instructions.
Expiration:	See kit insert for complete instructions.

Images




Western Blot

2D Western Blot of anti-E.coli HCP antibody. Load: 35 μ g
 Total HCP. pH: 5-8. Primary antibody: Rabbit anti-HCP
 antibody at 1:200 for overnight at 4°C. Secondary antibody:
 Goat anti-rabbit secondary antibody at 1:10,000 for 30 min
 at RT. Block: MB-070 for 1 hour at RT.



Kit Box

This product is assembled as a kit. See attached protocol or
 CofA for further details.



Western Blot

Western Blot of Rabbit anti-HCP antibody. Lane 1:
 Molecular Weight. Lane 2: Total HCP. Load: 10 μ g per
 lane. Primary antibody: Rabbit anti-HCP cocktail at 1:50 for
 overnight at 4°C. Secondary antibody: Goat anti-rabbit
 secondary antibody at 1:10,000 for 30 min at RT. Block: MB-
 070 for 1 hour at RT.

Disclaimer

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 5199, Limerick, Pennsylvania, USA.