

Datasheet for KHH001

Polymount Mounting Media MaxTag™ Histo for IHC**Overview**

Description:	Polymount Mounting Media MaxTag™ Histo for IHC - KHH001
Item No.:	KHH001
Size:	15 mL
Applications:	IHC

Product Details

Background:	PolyMount Mounting Media is used to prepare microscope slides and preserve samples under cover glass. PolyMount Mounting Media is ideal for investigators in Microbiology, Cancer, Immunology, Cell Biology, and Molecular Biology research.
Synonyms:	Polymount Mounting Media Maxtag™ Histo For IHC, Polymount Mounting Media For Immunohistochemistry
Detection Kit Type:	Immunohistochemistry Kit

Target Details

Relevant Links:	KHH001 SDS
------------------------	--

Application Details

Tested Applications:	IHC
Application Note:	This component is ready-to-use and requires no additional preparation.
Assay Dilutions:	All assays should be optimized by the user. Recommended dilutions (if any) may be listed below.
IF:	See kit insert.
IHC:	See kit insert.

Formulation

Physical State:	Liquid
Concentration:	Neat
Preservative:	None
Stabilizer:	None

Shipping & Handling

Shipping Condition:	Ambient
Storage Condition:	See kit insert for complete instructions.
Expiration:	See kit insert for complete instructions.

Disclaimer

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 5199, Limerick, Pennsylvania, USA.

