

Datasheet for KOA0103**Ms-Rb ALP SABC IHC Detection Kit****Overview**

Description:	Mouse and Rabbit Specific Alk Phos SABC Detection Kit - KOA0103
Item No.:	KOA0103
Size:	1 Kit
Reactivity:	Leporidae (Rabbit), Mus musculus (Mouse)
Host Species:	Goat

Product Details

Background:	Mouse and Rabbit IgG Specific Alkaline Phosphate Streptavidin-Biotin Complex Immunohistochemistry Detection Kit is a streptavidin-biotin immuno-enzymatic antigen detection system. This system involves the sequential incubation of the sample with an unconjugated primary antibody specific to the target antigen, a biotinylated secondary antibody that reacts with the primary antibody, enzyme-labeled streptavidin, and chromogen substrate. This kit includes reagents specifically designed for the use in the analysis of antigens in IHC and other immunodetection assays. The provided protocol includes steps for immunohistochemical analysis with low background using 5% BSA blocking reagent, biotinylated secondary antibody, an enzyme-labeled streptavidin, and a chromogenic agent.
Synonyms:	Alkaline Phosphatase/3,3'-diaminobenzidine Kit, StreptAvidin-Biotin Complex, biotin assay, Anti-Mouse IgG antibody, Anti-Rabbit IgG antibody, goat anti-Mouse/Rabbit IgG antibody, colormetric staining, immunohistochemistry staining
Host Species:	Goat
Conjugate:	Alkaline Phosphatase (AP)
Clonality:	Polyclonal
Detection Kit Type:	Immunohistochemistry Kit

Target Details

Reactivity:	Leporidae (Rabbit), Mus musculus (Mouse)
Relevant Links:	KOA0103 protocol

Application Details

Application Note:	This kit is ready to use and no dilution is needed. See kit insert for details. Useful for Immunohistochemistry.
Assay Dilutions:	All assays should be optimized by the user. Recommended dilutions (if any) may be listed below.
IHC:	User Optimized
Other:	Expiration Date: 07APR2023

Formulation

Shipping & Handling

Shipping Condition:	Wet Ice
Storage Condition:	Store at 4°C for 12 months. Avoid freezing. See kit insert for complete instructions.
Expiration:	Expiration date is one (1) year from date of receipt.

Disclaimer

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 5199, Limerick, Pennsylvania, USA.

