

Datasheet for KOA0141

Human Activin A ELISA Kit**Overview**

| | |
|----------------------|--|
| Description: | Human Activin A AccuSignal ELISA Kit - KOA0141 |
| Item No.: | KOA0141 |
| Size: | 1 Kit |
| Applications: | ELISA |
| Reactivity: | Human |
| Range: | 15.6 pg/ml - 1000 pg/ml |

Product Details

Background: Activin A is a homodimer of 14kDa beta-A. Activin A, a cytokine member of the transforming growth factor-beta superfamily, is expressed locally by the mesenchymal component of the hemopoietic microenvironment. Its expression is regulated on the mRNA level by different cytokines, and the biological activity of the protein is tightly controlled by several inhibitory molecules. Inhibins and activins are members of the transforming growth factor beta superfamily and are known to modulate the growth and differentiation of several cell types. Inhibins and activins inhibit and activate, respectively, the secretion of follitropin by the pituitary gland. Inhibins/activins are involved in regulating a number of diverse functions such as hypothalamic and pituitary hormone secretion, gonadal hormone secretion, germ cell development and maturation, erythroid differentiation, insulin secretion, nerve cell survival, embryonic axial development or bone growth, depending on their subunit composition. Inhibins appear to oppose the functions of activins. The standard product used in this kit is recombinant Activin A, which is composed of two single chains of 116 amino acids with the molecular mass of 26kDa.

| | |
|----------------------------|--|
| Synonyms: | Activin beta A chain, erythroid differentiation factor, erythroid differentiation protein, follicle stimulating hormone releasing protein, FSH releasing |
| Detection Kit Type: | ELISA Kit |
| Detection Range: | 15.6 pg/ml - 1000 pg/ml |
| Sensitivity: | <12 pg/ml |

Target Details

Gene Name: INHBA

| | |
|----------------------------|--|
| Reactivity: | Human |
| Immunogen: | Expression system for standard: CHO; Immunogen sequence: G311-S426 |
| Purity/Specificity: | Natural and recombinant human Activin A. There is no detectable cross-reactivity with other relevant proteins. |
| Relevant Links: | <ul style="list-style-type: none">• KOA0141 Protocol• UniProtKB - A4D1W7• NCBI - NP_002183.1• GeneID - 3624 |

Application Details

| | |
|-----------------------------|--|
| Tested Applications: | ELISA |
| Application Note: | Useful in Sandwich ELISA for Quantitative Detection of Antigen. Aliquot 0.1ml per well of the 1000pg/ml, 500pg/ml, 250pg/ml, 125pg/ml, 62.5pg/ml, 31.2pg/ml, 15.6pg/ml human Activin A standard solutions into the precoated 96-well plate. Add 0.1ml of the sample diluent buffer into the control well (Zero well). Add 0.1ml of each properly diluted sample of human cell culture supernates, serum, plasma(heparin, EDTA) or saliva to each empty well. It is recommended that each human Activin A standard solution and each sample be measured in duplicate. |
| Assay Dilutions: | All assays should be optimized by the user. Recommended dilutions (if any) may be listed below. |

Formulation

| | |
|-----------------------|----------------|
| Anticoagulant: | Heparin Sodium |
|-----------------------|----------------|


Shipping & Handling

| | |
|----------------------------|---|
| Shipping Condition: | Wet Ice |
| Storage Condition: | Store vials at 4°C prior to opening. Centrifuge product if not completely clear after standing at room temperature. This product is stable for 6 months at 4°C as an undiluted liquid. Dilute only prior to immediate use. For extended storage freeze at -20°C or below for 12 months. Avoid cycles of freezing and thawing. |
| Expiration: | See kit insert for complete instructions. |

Images

ELISA

Human Activin A AccuSignal ELISA Kit standard curve. Assay Range: 15.6pg/ml-1000pg/ml. Sensitivity: <12pg/ml.

**Disclaimer**

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 5199, Limerick, Pennsylvania, USA.