

Datasheet for KOA0205**Human IGFBP4 ELISA Kit****Overview**

Description:	Human IGFBP4 AccuSignal ELISA Kit - KOA0205
Item No.:	KOA0205
Size:	1 Kit
Applications:	ELISA
Reactivity:	Human
Range:	500 pg/ml - 32000 pg/ml

Product Details

Background:	Insulin-like growth factor-binding protein 4 is a protein that in humans is encoded by the IGFBP4 gene. This gene is a member of the insulin-like growth factor binding protein (IGFBP) family and encodes a protein with an IGFBP domain and a thyroglobulin type-I domain. The protein binds both insulin-like growth factors (IGFs) I and II and circulates in the plasma in both glycosylated and non-glycosylated forms. Binding of this protein prolongs the half-life of the IGFs and alters their interaction with cell surface receptors. In addition, IGFBP-4 is a unique protein and it consistently inhibits several cancer cells in vivo and in vitro. Its inhibitory action has been shown in vivo in prostate and colon. It is secreted by all colon cancer cells.
Synonyms:	AI875747, BP 4, BP4, Deb2, HT29 IGFBP, IBP 4, IBP-4, IBP4, IBP4_HUMAN, IGF binding protein 4, IGF-binding protein 4, IGF-BP4, IGFBP 4, IGFBP-4, IGFBP4, Insulin Like Growth Factor Binding Protein 4, Insulin-like growth factor-binding protein 4
Detection Kit Type:	ELISA Kit
Detection Range:	500 pg/ml - 32000 pg/ml
Sensitivity:	<20 pg/ml

Target Details

Gene Name:	IGFBP4
Reactivity:	Human
Immunogen:	Expression system for standard: NSO; Immunogen sequence: D22-E258

Purity/Specificity: Natural and recombinant human IGFBP4. There is no detectable cross-reactivity with other relevant proteins.

Relevant Links:

- [UniProtKB - P22692](#)
- [NCBI - NP_001543.2](#)
- [GeneID - 3487](#)
- [KOA0205 Protocol](#)

Application Details

Tested Applications: ELISA

Application Note: Useful in Sandwich ELISA for Quantitative Detection of Antigen. Aliquot 0.1ml per well of the 32000pg/ml, 16000pg/ml, 8000pg/ml, 4000pg/ml, 2000pg/ml, 1000pg/ml, 500pg/ml human IGFBP4 standard solutions into the precoated 96-well plate. Add 0.1ml of the sample diluent buffer into the control well (Zero well). Add 0.1ml of each properly diluted sample of human cell culture supernates, serum or plasma (heparin, EDTA) to each empty well. It is recommended that each human IGFBP4 standard solution and each sample be measured in duplicate.

Assay Dilutions: All assays should be optimized by the user. Recommended dilutions (if any) may be listed below.

Formulation

Anticoagulant: Heparin Sodium

Shipping & Handling

Shipping Condition: Wet Ice


Storage Condition: Store vials at 4°C prior to opening. Centrifuge product if not completely clear after standing at room temperature. This product is stable for 6 months at 4°C as an undiluted liquid. Dilute only prior to immediate use. For extended storage freeze at -20°C or below for 12 months. Avoid cycles of freezing and thawing.

Expiration: See kit insert for complete instructions.

Images

ELISA

Human IGFBP4 AccuSignal ELISA Kit standard curve. Assay
Range: 500pg/ml-32000pg/ml. Sensitivity: <20pg/ml.

**Disclaimer**

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 5199, Limerick, Pennsylvania, USA.