

**Datasheet for KOA0478****Human CSF1R - M-CSFR ELISA Kit****Overview**

|                      |                                                   |
|----------------------|---------------------------------------------------|
| <b>Description:</b>  | Human CSF1R/M-CSFR AccuSignal ELISA Kit - KOA0478 |
| <b>Item No.:</b>     | KOA0478                                           |
| <b>Size:</b>         | 1 Kit                                             |
| <b>Applications:</b> | ELISA                                             |
| <b>Reactivity:</b>   | Human                                             |
| <b>Range:</b>        | 62.5 pg/ml - 4000 pg/ml                           |

**Product Details**

|                            |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
|----------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Background:</b>         | CSF1R(Colony stimulating factor 1 receptor),also known as M-CSFR and CD115, is a cell-surface protein encoded, in humans, by the CSF1R gene. The gene is located on long arm of chromosome 5(5q32) on the Crick(minus) strand. The encoded protein is a tyrosine kinase transmembrane receptor and member of the CSF1/PDGF receptor family of tyrosine-protein kinases. The encoded protein is a single pass type I membrane protein and acts as the receptor for colony stimulating factor 1, a cytokine which controls the production, differentiation, and function of macrophages. Both CSF1R, and its ligand colony stimulating factor 1 play an important role in the development of the mammary gland and may be involved in the process of mammary gland carcinogenesis. |
| <b>Synonyms:</b>           | C FMS, CD 115, CD115, CD115 antigen, CFMS, Colony stimulating factor 1 receptor, CSF-1 receptor, CSF-1-R, CSF1R_HUMAN, CSFR, FIM2, FMS, FMS proto oncogene, FMS protooncogene, HDLS, M-CSF Receptor, M-CSF-R, Macrophage colony stimulating factor 1 receptor, Macrophage colony-stimulating factor 1 receptor, McDonough feline sarcoma viral (v fms) oncogene homolog, Oncogen FMS, Proto-oncogene c-Fms, V-FMS McDonough feline sarcoma viral oncogen homolog                                                                                                                                                                                                                                                                                                                 |
| <b>Detection Kit Type:</b> | ELISA Kit                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| <b>Detection Range:</b>    | 62.5 pg/ml - 4000 pg/ml                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| <b>Sensitivity:</b>        | <10 pg/ml                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |

**Target Details**

|                   |       |
|-------------------|-------|
| <b>Gene Name:</b> | CSF1R |
|-------------------|-------|

|                            |                                                                                                                                                                                                                         |
|----------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Reactivity:</b>         | Human                                                                                                                                                                                                                   |
| <b>Immunogen:</b>          | Expression system for standard: NSO; Immunogen sequence: I20-E512                                                                                                                                                       |
| <b>Purity/Specificity:</b> | Natural and recombinant human CSF1R. There is no detectable cross-reactivity with other relevant proteins.                                                                                                              |
| <b>Relevant Links:</b>     | <ul style="list-style-type: none"><li>• <a href="#">UniProtKB - P07333</a></li><li>• <a href="#">NCBI - NP_001275634.1</a></li><li>• <a href="#">GenelD - 1436</a></li><li>• <a href="#">KOA0478 Protocol</a></li></ul> |

## Application Details

|                             |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
|-----------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Tested Applications:</b> | ELISA                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Application Note:</b>    | Useful in Sandwich ELISA for Quantitative Detection of Antigen. Aliquot 0.1ml per well of the 4000pg/ml, 2000pg/ml, 1000pg/ml, 500pg/ml, 250pg/ml, 125pg/ml, 62.5pg/ml human CSF1R standard solutions into the precoated 96-well plate. Add 0.1ml of the sample diluent buffer into the control well (Zero well). Add 0.1ml of each properly diluted sample of human cell culture supernates, serum or plasma (heparin, EDTA) to each empty well. It is recommended that each human CSF1R standard solution and each sample be measured in duplicate. |
| <b>Assay Dilutions:</b>     | All assays should be optimized by the user. Recommended dilutions (if any) may be listed below.                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>ELISA:</b>               | 62.5pg/ml-4000pg/ml                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>Other:</b>               | DOE:<br>2-8° C - 08/20/2024<br>-20° C - 02/22/2025                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |

## Formulation


|                       |                |
|-----------------------|----------------|
| <b>Concentration:</b> | 1 Kit          |
| <b>Anticoagulant:</b> | Heparin Sodium |

## Shipping & Handling

|                            |         |
|----------------------------|---------|
| <b>Shipping Condition:</b> | Wet Ice |
|----------------------------|---------|

|                           |                                                                                                                                                                                                                                                                                                                               |
|---------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Storage Condition:</b> | Store vials at 4°C prior to opening. Centrifuge product if not completely clear after standing at room temperature. This product is stable for 6 months at 4°C as an undiluted liquid. Dilute only prior to immediate use. For extended storage freeze at -20°C or below for 12 months. Avoid cycles of freezing and thawing. |
| <b>Expiration:</b>        | See kit insert for complete instructions.                                                                                                                                                                                                                                                                                     |

## Images



### ELISA

Human CSF1R/M-CSFR AccuSignal ELISA Kit standard curve.  
Assay Range: 62.5pg/ml-4000pg/ml. Sensitivity: <10pg/ml.

## Disclaimer

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 5199, Limerick, Pennsylvania, USA.