

Datasheet for KOA0716**Mouse CD56 - NCAM1 ELISA Kit****Overview**

Description:	Mouse CD56/NCAM1 AccuSignal ELISA Kit - KOA0716
Item No.:	KOA0716
Size:	1 Kit
Applications:	ELISA
Reactivity:	Mouse
Range:	625 pg/ml - 40,000 pg/ml

Product Details

Background:	Neural Cell Adhesion Molecule, NCAM, also known as the cluster of differentiation CD56, is a hemophilic binding glycoprotein. It is a glycoprotein of Immunoglobulin(Ig) super family. At least 27 alternatively spliced NCAM mRNAs are produced, giving a wide diversity of NCAM isoforms. NCAM gene is located at 11q22-q23. This glycoprotein is mainly expressed on the surface of neurons, glia, skeletal muscle and natural killer cells. NCAM has been implicated as having a role in cell-cell adhesion, neurite outgrowth, synaptic plasticity, and learning and memory.
Synonyms:	CD56, N-CAM-1, NCAM-1, NCAM1, NCAM1_HUMAN, Neural cell adhesion molecule 1
Detection Kit Type:	ELISA Kit
Detection Range:	625 pg/ml - 40,000 pg/ml
Sensitivity:	<20 pg/ml

Target Details

Gene Name:	NCAM1
Reactivity:	Mouse
Immunogen:	Expression system for standard: NSO; Immunogen sequence: L20-T711
Purity/Specificity:	Natural and recombinant mouse CD56. There is no detectable cross-reactivity with other relevant proteins.
Relevant Links:	KOA0716 SDS UniProtKB - P13595

- [NCBI - NP_001074914.1](#)
- [GeneID - 17967](#)
- [KOA0716 Protocol](#)

Application Details

Tested Applications:	ELISA
Application Note:	Useful in Sandwich ELISA for Quantitative Detection of Antigen. Aliquot 0.1ml per well of the 40,000pg/ml, 20,000pg/ml, 10,000pg/ml, 5,000pg/ml, 2,500pg/ml, 1,250pg/ml, 625pg/ml mouse CD56 standard solutions into the precoated 96-well plate. Add 0.1ml of the sample diluent buffer into the control well (Zero well). Add 0.1ml of each properly diluted sample of mouse cell culture supernates, serum or plasma(heparin, EDTA) to each empty well. It is recommended that each mouse CD56 standard solution and each sample be measured in duplicate.
Assay Dilutions:	All assays should be optimized by the user. Recommended dilutions (if any) may be listed below.
ELISA:	625pg/ml-40,000pg/ml
Other:	DOE: 2-8° C - 08/07/2025 -20° C - 02/05/2026


Formulation

Concentration:	1 Kit
Anticoagulant:	Heparin Sodium

Shipping & Handling

Shipping Condition:	Wet Ice
Storage Condition:	Store vials at 4°C prior to opening. Centrifuge product if not completely clear after standing at room temperature. This product is stable for 6 months at 4°C as an undiluted liquid. Dilute only prior to immediate use. For extended storage freeze at -20°C or below for 12 months. Avoid cycles of freezing and thawing.
Expiration:	See kit insert for complete instructions.

Images



ELISA

Mouse CD56/NCAM1 ELISA Kit standard curve.

References

- Kagiava A et al. AAV9-mediated Schwann cell-targeted gene therapy rescues a model of demyelinating neuropathy. *Gene Ther.* (2021)

Disclaimer

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 5199, Limerick, Pennsylvania, USA.