

Datasheet for KOA0760**Human CRISP3 ELISA Kit****Overview**

Description:	Human CRISP3 AccuSignal ELISA Kit - KOA - KOA0760
Item No.:	KOA0760
Size:	1 Kit
Applications:	ELISA
Reactivity:	Human
Range:	1.56 ng/ml - 100 ng/ml

Product Details

Background:	Cysteine-rich secretory protein 3 is a cysteine-rich secretory protein that in humans is encoded by the CRISP3 gene. CRISP3 expression is found in the secretory epithelium throughout the male genital tract, with particularly high expression in the cauda epididymis and ampulla vas deferens. It is a potential biomarker for prostate cancer. CRISP3 gene can be defined as an early-transcribed gene with a putative role in the pathophysiology of sjogren's syndrome. CRISP3 is present in exocrine secretions and in secretory granules of neutrophilic granulocytes and is believed to play a role in innate immunity.
Synonyms:	Aeg 2, Aeg2, CRIS3_HUMAN, CRISP 3, CRISP-3, Crisp3, CRS 3, CRS3, Cysteine rich secretory protein 3, Cysteine-rich secretory protein 3, dJ442L6.3, MGC126588, OTTHUMP00000016590, OTTHUMP00000039911, SGP 28, SGP28, SGP28 protein, Specific granule protein (28 kDa), Specific granule protein of 28 kDa
Detection Kit Type:	ELISA Kit
Detection Range:	1.56 ng/ml - 100 ng/ml
Sensitivity:	<50 pg/ml

Target Details

Gene Name:	CRISP3
Reactivity:	Human
Immunogen:	Expression system for standard: NSO; Immunogen sequence: N21-Y245

Purity/Specificity: Natural and recombinant human CRISP3. There is no detectable cross-reactivity with other relevant proteins.

Relevant Links:

- [UniProtKB - P54108](#)
- [NCBI - NP_001177915.1](#)
- [GeneID - 10321](#)
- [KOA0760 Protocol](#)

Application Details

Tested Applications: ELISA

Application Note: Useful in Sandwich ELISA for Quantitative Detection of Antigen. Aliquot 0.1ml per well of the 100ng/ml, 50ng/ml, 25ng/ml, 1.25ng/ml, 6.25ng/ml, 3.12ng/ml, 1.56ng/ml human CRISP3 standard solutions into the precoated 96-well plate. Add 0.1ml of the sample diluent buffer into the control well (Zero well). Add 0.1ml of each properly diluted sample of human cell culture supernates, serum, plasma (heparin, EDTA), saliva, urine or tissue homogenates to each empty well. We recommend that each human CRISP3 standard solution and each sample is measured in duplicate.

Assay Dilutions: All assays should be optimized by the user. Recommended dilutions (if any) may be listed below.

Formulation

Anticoagulant: Heparin Sodium

Shipping & Handling

Shipping Condition: Wet Ice


Storage Condition: Store vials at 4°C prior to opening. Centrifuge product if not completely clear after standing at room temperature. This product is stable for 6 months at 4°C as an undiluted liquid. Dilute only prior to immediate use. For extended storage freeze at -20°C or below for 12 months. Avoid cycles of freezing and thawing.

Expiration: See kit insert for complete instructions.

Images

ELISA

Mouse Activin A AccuSignal Elisa Kit standard curve. Assay
Range: 15.6pg/ml-1000pg/ml. Sensitivity: <10pg/ml.

**Disclaimer**

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 5199, Limerick, Pennsylvania, USA.