

**Datasheet for KOA0763****Mouse CD6 ELISA Kit****Overview**

|                      |                                          |
|----------------------|------------------------------------------|
| <b>Description:</b>  | Mouse CD6 AccuSignal ELISA Kit - KOA0763 |
| <b>Item No.:</b>     | KOA0763                                  |
| <b>Size:</b>         | 1 Kit                                    |
| <b>Applications:</b> | ELISA                                    |
| <b>Reactivity:</b>   | Mouse                                    |
| <b>Range:</b>        | 62.5 pg/ml - 4000 pg/ml                  |

**Product Details**

|                            |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
|----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Background:</b>         | CD6(Cluster of Differentiation 6) is a human protein encoded by the CD6 gene. It is mapped to 11q12.2. CD6 is expressed predominantly in peripheral T cells and mature medullary thymocytes. This gene encodes a protein found on the outer membrane of T-lymphocytes as well as some other immune cells. The encoded protein contains three scavenger receptor cysteine-rich (SRCR) domains and a binding site for an activated leukocyte cell adhesion molecule. The gene product is important for continuation of T cell activation. What's more, this gene may be associated with susceptibility to multiple sclerosis. |
| <b>Synonyms:</b>           | CD_antigen=CD6, CD6, CD6 antigen Tp120, CD6 molecule, CD6_HUMAN, FLJ44171, OX52, T cell differentiation antigen CD6, T-cell differentiation antigen CD6, T12, TP120                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Detection Kit Type:</b> | ELISA Kit                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>Detection Range:</b>    | 62.5 pg/ml - 4000 pg/ml                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>Sensitivity:</b>        | <10 pg/ml                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |

**Target Details**

|                            |                                                                                                          |
|----------------------------|----------------------------------------------------------------------------------------------------------|
| <b>Gene Name:</b>          | CD6                                                                                                      |
| <b>Reactivity:</b>         | Mouse                                                                                                    |
| <b>Immunogen:</b>          | Expression system for standard: NSO; Immunogen sequence: L18-G396                                        |
| <b>Purity/Specificity:</b> | Natural and recombinant mouse CD6. There is no detectable cross-reactivity with other relevant proteins. |

**Relevant Links:**

- [UniProtKB - Q61003](#)
- [NCBI - NP\\_001032890.1](#)
- [GeneID - 12511](#)
- [KOA0763 Protocol](#)

## Application Details

|                             |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
|-----------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Tested Applications:</b> | ELISA                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Application Note:</b>    | Useful in Sandwich ELISA for Quantitative Detection of Antigen. Aliquot 0.1ml per well of the 4000pg/ml, 2000pg/ml, 1000pg/ml, 500pg/ml, 250pg/ml, 125pg/ml, 62.5pg/ml mouse CD6 standard solutions into the precoated 96-well plate. Add 0.1ml of the sample diluent buffer into the control well (Zero well). Add 0.1ml of each properly diluted sample of mouse cell culture supernates, serum or plasma(heparin, EDTA) to each empty well. It is recommended that each mouse CD6 standard solution and each sample be measured in duplicate. |
| <b>Assay Dilutions:</b>     | All assays should be optimized by the user. Recommended dilutions (if any) may be listed below.                                                                                                                                                                                                                                                                                                                                                                                                                                                  |

## Formulation

|                       |                |
|-----------------------|----------------|
| <b>Anticoagulant:</b> | Heparin Sodium |
|-----------------------|----------------|


## Shipping & Handling

|                            |                                                                                                                                                                                                                                                                                                                               |
|----------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Shipping Condition:</b> | Wet Ice                                                                                                                                                                                                                                                                                                                       |
| <b>Storage Condition:</b>  | Store vials at 4°C prior to opening. Centrifuge product if not completely clear after standing at room temperature. This product is stable for 6 months at 4°C as an undiluted liquid. Dilute only prior to immediate use. For extended storage freeze at -20°C or below for 12 months. Avoid cycles of freezing and thawing. |
| <b>Expiration:</b>         | See kit insert for complete instructions.                                                                                                                                                                                                                                                                                     |

## Images

**ELISA**

Mouse CD6 ELISA Kit standard curve.

**Disclaimer**

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 5199, Limerick, Pennsylvania, USA.