

Datasheet for MB-004**1.0 M Tris HCl pH 7.8****Overview**

Description:	1.0 M Tris HCl pH 7.8 - MB-004
Item No.:	MB-004
Size:	1 L

Product Details

Background:	Tris is useful in molecular biology and biochemistry due to its wide pH range (7-9) and its correspondence with the pH of living organisms. 1.0 M Tris HCl pH 7.8 is ideal for investigators in Molecular Biology, Cell Biology, and Immunology research.
Synonyms:	THAM, tris(hydroxymethyl)aminomethane

Target Details

Purity/Specificity:	1.0 M Tris HCl pH 7.8 was aseptically filtered through a Millipore 0.22 micron filter into clean, pre-sterilized containers. 1.0 M Tris HCl was tested on trypticase soy agar for 24 hours, 48 hours and 72 hours and was found to be negative for bacteria.
----------------------------	--

Application Details

Application Note:	1.0 M Tris HCl pH 7.8 is a concentrated stock solution and should be diluted appropriately with distilled, deionized water or equivalent to its final working concentration. This buffer consists of 1.0 M Tris HCl at pH 7.8. Meticulously prepared using ultra pure reagents dissolved in highly polished pharmaceutical grade deionized water.
Assay Dilutions:	All assays should be optimized by the user. Recommended dilutions (if any) may be listed below.

Formulation

Physical State:	Liquid (sterile filtered)
------------------------	---------------------------

Concentration:	1.0M
Buffer:	See application note.
Preservative:	None
Stabilizer:	None

Shipping & Handling

Shipping Condition:	Ambient
Storage Condition:	Store container at room temperature (18° to 26° C) prior to opening. If desired, the solution may be stored at 4° C or less. Some salts may precipitate out of solution at lower temperature. Allow buffer to equilibrate to room temperature (18° to 26° C) to restore solubility of some salts.
Expiration:	Expiration date is six (6) months from date of receipt.

Images



Bottle
1.0 M Tris HCl pH 7.8

References

- Buchbender A et al. Improved library preparation with the new iCLIP2 protocol. *Methods*. (2020)

Disclaimer

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 5199, Limerick, Pennsylvania, USA.