

Datasheet for RT-T093

Rat Liver

Overview

Description:	Rat Liver - RT-T093
Item No.:	RT-T093
Size:	1 Each
Origin:	Rat

Product Details

Species of Origin:	Rat
--------------------	-----

Target Details

Relevant Links:	RT-T093 SDS
-----------------	---

Formulation

Physical State:	Tissue
Buffer:	None
Sterility:	Non-sterile
Preservative:	None
Stabilizer:	None

Shipping & Handling

Shipping Condition:	Dry Ice
Storage Condition:	Store tissue at -20° C or colder prior to use.
Expiration:	No expiration date is given for this product if properly stored.

Disclaimer

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 5199, Limerick, Pennsylvania, USA.

